User's Manual

SERVO DUO GUARD

MAQUET
GETINGE GROUP

CRITICAL CARE

TABLE OF CONTENTS

| | | | |
|-----------|---|--|-----|
| 1 | Servo Duo Guard Instructions | | 4 |
| 2 | تعليمات Servo Duo Guard | | 9 |
| 3 | Servo Duo Guard Инструкции | | 14 |
| 4 | Servo Duo Guard-vejledning | | 19 |
| 5 | Servo Duo Guard juhised | | 24 |
| 6 | Servo Duo Guardian käyttöohjeet | | 29 |
| 7 | Instructions concernant le Servo Duo Guard | | 34 |
| 8 | Οδηγίες Servo Duo Guard | | 39 |
| 9 | Istruzioni per Servo Duo Guard | | 44 |
| 10 | サーボデュオガードの説明 | | 49 |
| 11 | Servo Duo Guard 滤菌器说明书 | | 54 |
| 12 | „Servo Duo Guard“ instrukcijos | | 59 |
| 13 | Instructies voor Servo Duo Guard | | 64 |
| 14 | Bruksanvisning for Servo Duo Guard | | 69 |
| 15 | Instrukcje dotyczące filtru Servo Duo Guard | | 74 |
| 16 | Instruções do Servo Duo Guard | | 79 |
| 17 | Instrukçioni Servo Duo Guard | | 84 |
| 18 | Инструкции по использованию фильтра Servo Duo Guard | | 89 |
| 19 | Uputstvo za Servo Duo Guard | | 94 |
| 20 | Pokyny pre Servo Duo Guard | | 99 |
| 21 | Instrucciones para Servo Duo Guard | | 104 |
| 22 | Instruktioner till Servo Duo Guard | | 109 |
| 23 | คำแนะนำการใช้งาน Servo Duo Guard | | 114 |
| 24 | Návod k použití filtru Servo Duo Guard | | 119 |
| 25 | Servo Duo Guard Talimatları | | 124 |
| 26 | Servo Duo Guard Gebrauchsanweisung | | 129 |
| 27 | Servo Duo Guard használati utasítások | | 134 |
| 28 | Hướng Dẫn Sử Dụng Servo Duo Guard | | 139 |


1 SERVO DUO GUARD INSTRUCTIONS

1.1 SYMBOLS

| Symbol | Explanation |
|---|---|
|  | Use by date |
|  | Do not re-use. Single use only. |
|  | Consult Instructions for use |
|  | Keep away from sunlight |
|  | Do not use if packaging is damaged |
| REF | Order number |
| LOT | Number to identify the production batch |
|  | Manufacturer |
|  0123 | CE label—indicates compliance with the requirements of the Medical Device Directive 93/42/EEC |

1.2 GENERAL INFORMATION

The Servo Duo Guard is a highly efficient single use bacterial and viral filter for application in respiratory care and anesthesia. The filter provides filtration (bacterial and viral efficiency of 99.9999%) for reducing possible cross contamination between patient and equipment.

The unique dual filter design of the Servo Duo Guard minimizes the risk of sudden increases in expiratory resistance, which is a known problem when using filters during nebulization. This is made possible by combining a high quality HEPA filter for bacterial/viral protection (1) with an electrostatic filter (2) for collecting nebulization residuals from the patient's expired gases.


The Servo Duo Guard is bi-directional and intended to be used on both the expiratory and/or inspiratory limbs of breathing circuits. If nebulization is used, the Servo Duo Guard must be positioned in a specific way described under Instructions for Use.

The filter housing is transparent to allow visual inspection while in use.

1.3 INSTRUCTIONS FOR USE

Place the filter either on the expiratory or inspiratory limb of the ventilator (or anesthetic breathing system).

Note: For nebulization, position the filter on the expiratory limb with the arrows (3) on the filter pointing in the gas flow direction.

The filter is designed to work during 48 hours.

For nebulization, the filter has been tested for the following drugs under the following conditions:

- Using a known drug intended for nebulization (e.g. Berodual®, Acetylcystein, Combivent, Garamycin, Pulmicort*, Bromhexin, Tacholiquin, Tyloxapol, DNCG Stada).

- Using the AERONEB or Ultrasonic nebulizers provided by MAQUET.
- Nebulizing a maximum of 6 times 5 ml (4 ml Pulmicort) during 48 hours, where each nebulization period is ≤ 30 minutes with continuous use of single or dual-heated active humidification (100% RH) between nebulization.

WARNINGS!

- For single patient use only.
- Handle used filter as contaminated waste.
- If nebulization is performed outside of the test frame described on the previous page, monitor the patient airway pressure with extra care. Also check for particle collection of, for instance, blood.
- If active humidification is used, continuously check to ensure that possible condensation created in the expiratory limb does not affect ventilator performance.

CAUTIONS:

- To be replaced every 48 hours.
- Placement of the filter between the tracheal or tracheostomy tube and the breathing circuit will increase the dead space by 170 ml. Only use this filter on the inspiratory or expiratory limb (machine ends) of the breathing circuit.


1.4 AFTER USE

- Discard the filter. The filter should be disposed of in accordance with hospital procedures for contaminated waste.
- Avoid spillage of any liquid remaining in the filter. The liquid may be contaminated.

1.5 ORDER INFORMATION

Order number 66 71 775 represents a carton of 60 units.

Expiry date yyyy-mm (1) and lot number (2) is stated on each Servo Duo Guard as shown in illustration below.

1.6 TECHNICAL DATA

1.6.1 GENERAL

The device is intended for adult, pediatric and infant populations.

Technical data

| | |
|---|--|
| Dimensions | Length = 110 mm, maximum diameter = 78 mm |
| Weight | 55 g |
| Internal volume | 170 ml |
| Filter medium | Hydrophobic HEPA filter (High Efficiency Particulate Air) Electrostatic filter |
| Filter area | 955 cm ² |
| Flow range | 0-3.3 l/s |
| Pressure | Max. 30 kPa |
| Leakage | <2 ml/min at 15 kPa over pressure during 5 seconds |
| Bacterial (BFE) efficiency ¹ | 99.9999% |
| Viral (VFE) efficiency ² | 99.9999% |
| Quality check | ISO 23328-1:2003 Each filter is tested individually for filter effectiveness and leakage. |
| Manufacture | Hygienically manufactured and packed. |
| Environment | Incineration of the device will produce carbon dioxide and water. |

Technical data

| | |
|-----------------------------------|---|
| Storage and transport temperature | Room temperature. For short periods (1 week): Ambient temperature -25°C to +60°C. Relative humidity < 95%. |
| Shelf life | 5 years from manufacturing date, when properly stored in unbroken package |
| Material | Thermoplastic housing, hydrophobic glass fibre paper |

REFERENCES

¹Nelson Laboratories. Bacterial Filtration Efficiency Test BFE No 350457, with test values >99.9999%

²Nelson Laboratories. Virus Filtration Efficiency Test VFE No 350467, with test values >99.9999%

1.6.2 RESISTANCE

| | |
|------------|---|
| Dry filter | <75 Pa (0.75 cmH ₂ O) at 0.5 l/s |
| | <150 Pa (1.5 cmH ₂ O) at 1 l/s |
| | <330 Pa (3.3 cmH ₂ O) at 2 l/s |
| Wet filter | <100 Pa (1.0 cmH ₂ O) at 0.5 l/s |
| | <200 Pa (2.0 cmH ₂ O) at 1 l/s |
| | <450 Pa (4.5 cmH ₂ O) at 2 l/s |
| | (Wet filter = 24 hours use at 10 l/min and 100% RH) |

1.6.3 CONNECTIONS

| | |
|------------------------------|---------------------|
| Coaxial cone male/ female | ISO M22/F15 (inlet) |
| Cone female | ISO F22 (outlet) |

1.7 CERTIFICATES

1.7.1 CHINA

| | |
|--|--|
| Certificate number: | SFDA (I) 2012 2660032 |
| Registration product standard number: | YZB/SWE 4734-2011 (Filter) |
| After sales service agent: | MAQUET (Shanghai) Medical Equipment Co., Ltd. |
| After sales service address: | Suite D6 (01-21) & D7, Floor 1, No. 56 Meisheng Road, Waigaoqiao Free Trade Zone, Shanghai, China |
| After sales service telephone: | 800-820-0207 |

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For local contact:
Please visit our website
www.maquet.com

0123

GETINGE GROUP is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, GETINGE and MAQUET. ArjoHuntleigh focuses on patient mobility and wound management solutions. GETINGE provides solutions for infection control within healthcare and contamination prevention within life sciences. MAQUET specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.

SERVO DUO GUARD

User's Manual

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Maquet Critical Care AB