The Bard logo consists of the word "BARD" in a bold, black, stylized font. The letters are blocky with some internal cutouts, particularly in the 'A' and 'R's.The CE mark is a large, bold, black 'C' followed by a large, bold, black 'E'. Below the 'E' is the number '0086' in a smaller, bold, black font.The logo for Bard Access Systems, Inc. features a small icon of a building with a flag on top, followed by the text "Bard Access Systems, Inc." in a bold, black, sans-serif font.

605 North 5600 West
Salt Lake City, Utah 84116 U.S.A.
801-522-5000

Customer Service: 800-545-0890
Clinical Information: 800-443-3385

The website URL "www.bardaccess.com" is displayed in a bold, black, sans-serif font. It is enclosed within a stylized graphic element that resembles a curved arrow or a swoosh, with a red-to-black gradient.The EC REP logo consists of the letters "EC" in a bold, black, sans-serif font, followed by the letters "REP" in a bold, black, sans-serif font. Both are enclosed within a black rectangular border.

Bard Limited
Forest House, Brighton Road
Crawley, West Sussex
RH11 9BP UK

0728531
1104R

Hickman*, Leonard* and Broviac* Central Venous Catheters — Long Term

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











Rx Only

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

La loi fédérale (E.-U.) n'autorise la délivrance de ce système qu'à un médecin ou sur prescription d'un médecin.
Gemäß den US-Bundesgesetzen darf dieses Produkt nur von einem Arzt oder auf Anordnung eines Arztes verkauft werden.
In base alla legge federale (U.S.A.) La vendita di questo dispositivo è consentita solo da parte o su prescrizione di un medico.
La ley federal (EE.UU.) limita la venta de este dispositivo, que debe ser realizada por un médico o por prescripción facultativa.
Krachtens de federale wetgeving (VS) mag dit instrument uitsluitend door of op voorschrift van een arts worden verkocht.
As Leis Federais dos E.U. determinam que este dispositivo seja vendido somente por médicos ou por ordem destes.
Το Ομοσπονδιακό Δίκαιο των Η.Π.Α. περιορίζει την πώληση αυτής της συσκευής μόνο από γιατρούς ή κατόπιν εντολής γιατρού.
Federal lov i USA begrænser denne anordning til salg af eller på bestilling af en læge.
Enligt federal lagstiftning (USA) får denna produkt säljas endast av läkare eller enligt läkares ordination.
Hydysvaltain lain mukaan tätä tuotetta saa myydä vain lääkärit tai lääkärin määräyksestä.
Amerikansk lovgivning begrænser dette udstyret til salg af lege eller etter henvisning av lege.
Ezen eszköz kizárólag orvos által, illetve orvosi rendelvényre értékesíthető, az Egyesült Államok vonatkozó szövetségi törvénye értelmében.
Federální zákony USA omezují prodej tohoto zařízení jen na prodej lékařům nebo prodej na objednávku lékaře.
Zgodnie z prawem federalnym (w USA) urządzenie to może być sprzedawane wyłącznie przez lekarzy lub na ich zlecenie.
A.B.D. Federal Yasalarina göre bu cihaz bir hekim tarafından veya izniyle satılabilir.
Согласно федеральному закону (США), это устройство может быть продано только врачам или по их заказу.

**Open-Ended Catheter**

cathéter à bout ouvert
Endoffener Katheter
Catetere a punta aperta
Catéter de extremo abierto
Open-eind-catheter
Cateter aberto
Καθετήρας Ανοικτού Ακροῦ
Kateter med åben ende
Kateter med öppen spets
Avoim katetri
Kateter med åpen ende
Nýitott végú katéter
Katétr s otevřeným koncem
Cewnik otwarty
Açık uçlu katéter
Катетер с открытым концом

**French Size**

Taille en French
Größe (F)
Misura french
Tamaño French
French-formaat
Tamanho francês
Μέγεθος σε French
French-størrelse
Storlek i French
French-koko
French str.
Francia méret
Velikost (francouzská)
Rozmiar wg numeracji francuskiej
Boyutu (French cinsinden)
Французская шкала

**Triple Lumen**

à trois lumières
Dreilumig
Tripló lume
de tres lúmenes
Drievoudig lumen
Lúmen Tripló
Τριπλού Αιλού
Trippellumen
Trippellumen
Kolmitiehyeinen
Tredobbeltlumen
Háromüregű
Trojitý lumen
Trójdrożny
Üç lümen
Трёхпросветный катетер

**Single Lumen**

à une lumière
Einlumig
Singolo lume
de un solo lumen
Eén lumen
Lúmen Único
Μονού Αιλού
Enkeltlumen
Enkellumen
Yksitiehyeinen
Enkeltlumen
Együregű
Jednoduchý lumen
Jednodrożny
Tek lümen
Однопросветный катетер

**Peel-Apart**

pelable
Peel-Apart
Peel away
Pelable
Splijtbaar
Abra por Agua
Με Διαχωριζόμενα Τμήματα
Peel-apart
Isádragbar
Peel-apart
Splittes ("peel apart")
Széthúzható
Stáhnout
Rozdzierany
Ayrılabilir
Разделить

**Dual Lumen**

à deux lumières
Doppellumig
Doppio lume
de dos lúmenes
Dubbel lumen
Lúmen Duplo
Διπλού Αιλού
Dobbeltlumen
Dubbeltlumen
Kaksitiehyeinen
Dobbeltlumen
Kétüregű
Dvojity lumen
Dwudrożny
Çift lümen
Двухпросветный катетер

BAIRD**Access Systems**

Introduction

Description:


Hickman[®], Leonard[®], and Broviac[®] Central Venous Catheters are constructed of specially formulated and processed silicone. The catheters are radiopaque with female luer locking connectors and **SureCuff[®]** Tissue Ingrowth Cuffs for fixation of the catheters in a subcutaneous tunnel. Each catheter is provided in a double sterile package.

Placement:


The catheter is placed into one of the large central veins so the tip lies in the superior vena cava above the right atrium. It is tunneled subcutaneously to the desired exit site. The **SureCuff[®]** Tissue Ingrowth Cuff, attached to the catheter, is positioned 3-5 cms below the skin exit site in the tunnel. The cuff promotes tissue ingrowth to secure the catheter in place.

Schematics:


Single Lumen

Dual Lumen

Triple Lumen

Indications For Use

Hickman*, **Leonard*** and **Broviac*** Catheters are designed for long-term vascular access and for use in patients that lack adequate peripheral venous access. They are available in single, dual and triple lumen catheters.

All **Hickman***, **Leonard*** and **Broviac*** central venous catheters are designed for the administration of I.V. fluids, blood products, drugs, and parenteral nutrition solutions, as well as blood withdrawal.

Note: While smaller lumen **Broviac*** catheters have been used successfully for blood withdrawal, their small lumen sizes increase the chance of clotting.

VitaCuff* Antimicrobial Cuff


Description

The **VitaCuff*** device is designed to help provide protection against infections related to vascular access catheters. The outer, tissue-interfacing surface of the **VitaCuff*** device may help reduce the incidence of infection by incorporating an antimicrobial agent into the porous collagen matrix.

The **VitaCuff*** device is comprised of two concentric layers of material. The internal layer is constructed of specially formulated and processed medical grade silicone. The external, tissue-interfacing layer is **VitaGuard*** antimicrobial collagen matrix. The antimicrobial activity of the **VitaGuard*** material is attributable to the silver ions bound to the collagen matrix. The activity lasts until the **VitaGuard*** matrix is completely absorbed by the tissue in four to six weeks. (Refer to figure 1)

The **VitaGuard*** collagen sponge is initially in a compressed state for ease of insertion. After placement, the matrix absorbs physiological fluids, quickly expands to approximately twice its original size, and helps provide an antimicrobial barrier and a physical barrier at the exit site. Tissue ingrowth into the **VitaGuard*** collagen matrix occurs in a few days, further securing the catheter in place, and reducing catheter movement.

Proper VitaCuff* Positioning

Caution: The antimicrobial cuff is not intended to be used as a treatment for catheter related infections. The antimicrobial cuff does not provide protection against "blood seeding" infection or infusate-related infection. It is not intended to provide protection from bacteria for longer than one month. The antimicrobial cuff should not be used on patients with known sensitivities to silver ions or collagen.

Contraindications, Warnings, Cautions and Precautions

Contraindications


The device is contraindicated whenever:

- The presence of device related infection, bacteremia, or septicemia is known or suspected.
- The patient's body size is insufficient to accommodate the size of the implanted device.
- The patient is known or is suspected to be allergic to materials contained in the device.
- Severe chronic obstructive lung disease exists (percutaneous subclavian placement only).
- Past irradiation of prospective insertion site.
- Previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site.
- Local tissue factors will prevent proper device stabilization and/or access.
- Do not use the antimicrobial cuff in patients with known sensitivities to silver or collagen.

Warnings:

- Intended for Single Use. DO NOT REUSE. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.
- This is not a right atrium catheter. Avoid positioning the catheter tip in the right atrium. Placement or migration of the catheter tip into the right atrium may cause cardiac arrhythmia, myocardial erosion or cardiac tamponade. The risk of these potential complications may be more likely in neonatal patients.
- Avoid vessel perforation.
- Hold thumb over exposed orifice of sheath to prevent air aspiration. The risk of air aspiration is reduced by performing this part of the procedure with the patient performing the Valsalva maneuver.
- You should not feel any resistance when withdrawing the catheter from the vein. If you do encounter resistance, this may indicate that the catheter is being pinched between the clavicle and first rib (the "pinch-off" sign). Do not continue pulling against resistance as this may cause catheter breakage and embolism. Free up the resistance (e.g. by repositioning the patient) before proceeding further.
-  After use, this product may be a potential biohazard. Handle and discard in accordance with accepted medical practice and applicable local, state and federal laws and regulations.
- If the artery is entered, withdraw the needle and apply manual pressure for several minutes. If the pleural space is entered, withdraw the needle and evaluate patient for possible pneumothorax.

- **Pinch-off Prevention:** Catheters placed percutaneously or through a cut-down, into the subclavian vein, should be inserted at the junction of the outer and middle thirds of the clavicle, lateral to the thoracic outlet. The catheter should not be inserted into the subclavian vein medially, because such placement can lead to compression of the catheter between the first rib and the clavicle, which can cause damage and even severance of the catheter. A radiographic confirmation of catheter placement should be made to ensure that the catheter is not being pinched by the first rib and clavicle.^{1,2}

Signs of Pinch-off

Clinical:

- Difficulty with blood withdrawal
- Resistance to infusion of fluids
- Patient position changes required for infusion of fluids or blood withdrawal

Radiologic:

- Grade 1 or 2 distortion on chest X-ray.

Pinch-off should be evaluated for degree of severity prior to explantation. Patients indicating any degree of catheter distortion at the clavicle/first rib area should be followed diligently. There are grades of pinch-off that should be recognized with appropriate chest x-ray as follows:^{3,4}

Grade	Severity	Recommended Action
Grade 0	No distortion	No action.
Grade 1	Distortion present without luminal narrowing	Chest x-ray should be taken every one to three months to monitor progression of pinch off to grade 2 distortion. Shoulder positioning during chest x-rays should be noted as it can contribute to changes in distortion grades.
Grade 2	Distortion present with luminal narrowing	Removal of the catheter should be considered.
Grade 3	Catheter transection or fracture	Prompt removal of the catheter.

Cautions:

- Carefully read and follow all instructions prior to use.
- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Only qualified healthcare practitioners should insert, manipulate and remove these devices.
- When tunneling, the catheter must not be forced.
- Avoid inadvertent puncture of the skin or fascia with the tip of the tunneler.
- The entire collagen (tan) portion of the **VitaCuff*** Antimicrobial Cuff must be placed beneath the skin level to avoid migration of the cuff out of the tunnel and exit site.
- Do not insert guidewire beyond the bevel of the needle while removing straightener from the needle hub in order to prevent guidewire damage or shearing.
- If the guidewire must be withdrawn while the needle is inserted, remove both the needle and guidewire as a unit to help prevent the needle from damaging or shearing the guidewire.
- Do not grasp the catheter with any instrument that might sever or damage the catheter.
- Do not cut the catheter before removal from vein to avoid catheter embolism.
- Do not use scissors or any sharp-edged instruments as they could damage the catheter.

Precautions:

- Follow Universal Precautions when inserting and maintaining the catheter.
- Follow all contraindications, warnings, cautions, precautions and instructions for all infusates as specified by its manufacturer.
- Use aseptic techniques whenever the catheter lumen is opened or connected to other devices. Povidone-iodine is the suggested antiseptic to use with this device and components. Acetone and tincture of iodine should not be used because they could adversely affect the performance of the catheter and connectors. 10 % acetone / 70 % isopropyl alcohol swabsticks used for dressing changes should not adversely affect the catheter.

I. Prior to beginning placement procedure, do the following:

- Examine package carefully before opening to confirm its integrity and that the expiration date has not passed. The device is supplied in a double sterile package and is non-pyrogenic. Do not use if package is damaged, opened or the expiration date has passed. Sterilized by ethylene oxide. Do not Resterilize.
- Inspect kit for inclusion of all components.
- When device includes an antimicrobial cuff, do not expose the cuff to fluids prior to insertion. Handle carefully to avoid cuff damage.
- Fill (prime) the device with sterile heparinized saline or normal saline solution to help avoid air embolism.
- When using an introducer kit, verify that the catheter fits easily through the introducer sheath.

II. To avert device damage and/or patient injury during placement:

- Avoid accidental device contact with sharp instruments and mechanical damage to the catheter material. Use only smooth-edged atraumatic clamps or forceps.
- Avoid perforating, tearing or fracturing the catheter when using a guidewire.

- Do not use the catheter if there is any evidence of mechanical damage or leaking.
- Avoid sharp or acute angles during implantation which could compromise the patency of the catheter lumen(s).
- If sutures are used to secure the catheter, make sure they do not occlude or cut the catheter.
- When using percutaneous introducers:
 - Carefully insert the introducer and catheter to avoid inadvertent penetration to vital structures in the thorax.
 - To avoid blood vessel damage, do not allow the percutaneous introducer sheath to remain indwelling in the blood vessel without the internal support of a catheter or dilator.
 - Simultaneously advance the sheath and dilator with rotational motion to help prevent sheath damage.
- During insertion of catheter with antimicrobial cuff:
 - Minimize the exposure of the cuff to pooled blood by sponging the intended cuff placement site.
 - The entire collagen (tan) portion of the cuff **must** be in the subcutaneous tissue at the catheter exit site.

III. After placement, observe the following precautions to avoid device damage and/or patient injury:

- Do not use the catheter if there is any evidence of mechanical damage or leaking. Damage to the catheter may lead to rupture, fragmentation and possible embolism and surgical removal.
- Accessories and components used in conjunction with this device should incorporate Luer lock connections.
- If signs of extravasation exist, discontinue injections. Begin appropriate medical intervention immediately.
- Infusion pressure greater than 25 psi (172 kPa) may damage blood vessels and viscous and is not recommended. DO NOT USE A SYRINGE SMALLER THAN 10 ml!

Possible Complications

The use of an indwelling central venous catheter provides an important means of venous access for critically ill patients; however, the potential exists for serious complications including the following:

- Air Embolism
- Allergic Reaction to Silver or Collagen (Catheters with **VitaCuff*** Antimicrobial Cuff only)
- Bleeding
- Brachial Plexus Injury
- Cardiac Arrhythmia
- Cardiac Tamponade
- Catheter or Cuff Erosion Through Skin
- Catheter Embolism
- Catheter or Cuff Occlusion
- Catheter Occlusion, Damage or Breakage due to Compression Between the Clavicle and First Rib
- Catheter-related Sepsis
- Endocarditis
- Exit Site Infection
- Exit Site Necrosis
- Extravasation
- Fibrin Sheath Formation
- Hematoma
- Hemothorax
- Hydrothorax
- Intolerance Reaction to Implanted Device
- Laceration of Vessels or Viscus
- Myocardial Erosion
- Perforation of Vessels or Viscus
- Pneumothorax
- Spontaneous Catheter Tip Malposition or Retraction
- Thoracic Duct Injury
- Thromboembolism
- Venous Thrombosis
- Ventricular Thrombosis
- Vessel Erosion
- Risks Normally Associated with Local and General Anesthesia, Surgery, and Post-Operative Recovery

These and other complications are well documented in medical literature and should be carefully considered before placing the catheter.

Hickman*, Leonard* and Broviac* Central Venous Catheter Placement Procedures

Section A: Prepping Procedure

Before beginning procedure, read the “Contraindications, Warnings, Cautions and Precautions” and “Possible Complications” sections of this manual.

1. Create sterile field and open tray.
2. Prep venipuncture/cutdown area, tunnel and tunnel exit areas.
3. Perform local anesthetic infiltration in venipuncture/cutdown, tunnel and tunnel exit site areas.
4. Irrigate the catheter with sterile heparinized saline (100 u/ml) and inspect for leakage. Clamp the catheter over the clamping sleeve(s).
5. Place patient in the Trendelenburg position with head turned away from the intended venipuncture site.

Refer to Section C or D Prior to Tunneling Procedure

Section B: Tunneling Procedure

1. Measure catheter against chest wall of patient to determine desired location of **SureCuff*** Tissue Ingrowth Cuff and exit site. Mark locations.
2. Tunneling procedure.

Note: The subcutaneous tunnel should be approximately 10 to 15 cm long with the **SureCuff*** Tissue Ingrowth Cuff positioned in the tunnel. The cuff will be less prominent if positioned over an intercostal space.

Single Lumen Catheters

Create subcutaneous tunnel from venous entrance site to skin exit site using tunneler or long forceps. (Refer to figure 2)

- a. Advance the tip of the tunneler from the venous entry site down to the desired catheter exit site.
- b. Thread suture material through the suture eye and tie it around the catheter tip when the tip of the tunneler emerges through the exit site.
- c. Pull the catheter up through the tunnel to the venous entry site. (Initial resistance may be felt as the **SureCuff*** Tissue Ingrowth Cuff or **VitaCuff*** Antimicrobial cuff first enters the tunnel.) Gently holding the catheter proximal to the cuff while pulling the tunneler and catheter through the subcutaneous tunnel should result in smooth passage of the cuff into the tunnel. **Caution:** When tunneling, the catheter must not be forced.
- d. Cut off the end of the catheter tied by suture.

- e. Estimate the catheter length required for the tip placement at the junction of the superior vena cava and right atrium by placing the catheter on the chest along the venous path to the right atrium. Cut catheter to length at a 45° angle.

Multi-Lumen Catheters:


Create subcutaneous tunnel from skin exit site to venous entrance using tunneler or long forceps. (Refer to figure 3)

- a. Grasp the tunneler at the end with protective cover.
- b. Insert the rounded tip of the tunneler into a small incision at the desired catheter exit site.
- c. Form tunnel by advancing the tip of the tunneler from the skin exit site up to the venous entry site.
Caution: Avoid inadvertent puncture of the skin or fascia with the tip of the tunneler.
- d. Remove the protective cover and attach one of the lumen tips onto the tunneler barb with a twisting motion. Barb threads must be completely covered by the catheter tip to adequately secure the catheter as it is pulled through the tunnel. A suture may be tied around the catheter between the tunneler body and large barb to hold it more securely.
- e. Pull the catheter up through the tunnel to the venous entry site. (Initial resistance may be met as the **VitaCuff®** Antimicrobial Cuff or **SureCuff®** Tissue Ingrowth Cuff first enters the tunnel.) Gently holding the catheter proximal to the cuff while pulling the tunneler and catheter through the subcutaneous tunnel should result in smooth passage of the cuff into the subcutaneous tunnel. **Caution:** When tunneling, the catheter must not be forced.
Caution: The entire collagen (tan) portion of the **VitaCuff®** Antimicrobial Cuff must be placed beneath the skin level to avoid migration of the cuff out of the tunnel and exit site.
- f. Remove the catheter tip from the tunneler barb.
- g. Cut the catheter to length at a 45° angle with the smaller lumen cut shorter than the larger lumen.

For Percutaneous Placement see Section D.

Section C: Cutdown Technique

1. Surgically isolate the desired vessel through a small skin incision.
Note: The external jugular vein, the cephalic vein at the delto-pectoral groove, and the axillary subclavian vein are the most common vessels used for catheter insertion. It may be necessary to use the internal jugular vein for insertion of larger catheters.
2. Refer to section B for catheter measurement and tunneling procedure.
3. Insert the catheter through a small venotomy into the isolated vein and advance to desired position in vessel.
4. Verify catheter tip location radiographically. The preferred location of the catheter tip is at the junction of the superior vena cava and the right atrium.
Warning: This is not a right atrium catheter. Avoid positioning the catheter tip in the right atrium. Placement or migration of the catheter tip into the right atrium may cause cardiac arrhythmia, myocardial erosion or cardiac tamponade. The risk of these potential complications may be more likely in neonatal patients.

5. Unclamp catheter and draw blood through the lumen(s) of the catheter to insure patency after placement is complete, but before closing the skin at the venipuncture site. If catheter is not patent, adjust catheter at curvature point to relieve possible restriction. Irrigate catheter lumen(s) with 10 ml of normal saline to clear catheter of blood. Instill sterile heparinized saline per lumen to create a heparin lock. Clamp catheter.
6. Attach injection cap(s) or connect to intravenous fluid source.
7. Close the skin at the venipuncture site as necessary, taking care not to damage the catheter.
8. Suture catheter at exit site.
9. Secure catheter at exit site with a sterile dressing. The external segment of the catheter should be coiled and taped. Avoid tension on the catheter segment to prevent dislodging the catheter.

Section D: Percutaneous Technique

Before beginning procedure, read the “Contraindications, Warnings, Cautions and Precautions” and “Possible Complications” sections of this manual.

1. Locate desired vessel using a small needle attached to a syringe. **Note:** The subclavian vein is entered percutaneously at the point that identifies the junction of the outer and middle thirds of the clavicle using the needle and syringe. (Refer to figure 4)

Refer to the “Warnings” section concerning Catheter Pinch-off.

2. Attach introducer needle to the syringe and insert into vessel alongside the small needle. Remove small needle.
3. Aspirate gently as the insertion is made. **Warning:** If the artery is entered, withdraw the needle and apply manual pressure for several minutes. If the pleural space is entered, withdraw the needle and evaluate patient for possible pneumothorax. (Refer to figure 5)
4. When the subclavian vein has been entered, remove the syringe leaving the needle in place. Place a finger over the hub of the needle to minimize blood loss and the risk of air aspiration. The risk of air aspiration is reduced by performing this part of the procedure with the patient performing the Valsalva maneuver. (Refer to figure 6)
5. Straighten “J” tip of guidewire with tip straightener and insert tapered end of tip straightener into the needle. Tip straightener should not be advanced over the guidewire beyond the guidewire tip. **Caution:** Do not insert guidewire beyond the bevel of the needle while removing straightener from the needle hub in order to prevent guidewire damage or shearing. Remove the tip straightener and advance the guidewire into the superior vena cava. Advance the guidewire as far as appropriate for the procedure. Verify correct positioning radiographically. (Refer to figure 7)
6. Gently withdraw and remove needle. **Caution:** If the guidewire must be withdrawn while the needle is inserted, remove both the needle and guidewire as a unit to help prevent the needle from damaging or shearing the guidewire. (Refer to figure 8)

- Refer to section B for catheter measurement and tunneling procedure.
- Make a small (approx. 1 cm wide) incision parallel to the clavicle, positioning the guidewire at the center of the incision to permit proper entry of vessel dilator and sheath introducer.

Intro-Eze* Introducer Instruction

(For Peel-Apart Introducer see #16)

- Advance the vessel dilator and sheath introducer as a unit over the exposed guidewire using a rotational motion. Advance it into the subclavian vein as a unit, leaving at least 2 cms of sheath exposed. **Warning:** Avoid vessel perforation. (Refer to figure 9)
- Withdraw the vessel dilator and "J" guidewire, leaving the sheath in place. **Warning:** Hold thumb over exposed orifice of sheath to prevent air aspiration. The risk of air aspiration is reduced by performing this part of the procedure with the patient performing the Valsalva maneuver. (Refer to figure 10)
- Advance the catheter through the sheath and into the vein. (Refer to figure 11)
- Verify catheter tip location radiographically. **Warning:** This is not a right atrium catheter. Avoid positioning the catheter tip in the right atrium. Placement or migration of the catheter tip into the right atrium may cause cardiac arrhythmia, myocardial erosion or cardiac tamponade. The risk of these potential complications may be more likely in neonatal patients. Preferred location of the catheter tip is at the junction of the superior vena cava and the right atrium.
- Pull the storage tube from the slit. Place the channeled portion of the slit onto the catheter near the proximal end of the introducer sheath.
- Grasp the proximal end of the slit between the thumb and index finger of one hand. With the tips of the fingers, reach around the slit and secure the catheter into the channeled portion. (Refer to figure 12)
- Withdraw the sheath over the catheter, sliding the proximal opening of the sheath over the nose of the channel and into the blade. Continue to withdraw the sheath, pulling it away from the catheter, until it is completely slit. Remove and discard the slit sheath and slit. (Refer to figure 13)

Proceed to step 23

Peel-Apart Sheath Introducer Instructions:

- Advance the vessel dilator and sheath introducer as a unit over the exposed guidewire using a rotational motion. Advance it into the subclavian vein as a unit, leaving at least 2 cms of sheath exposed. **Warning:** Avoid vessel perforation. (Refer to figure 14)
- Squeeze the hub handles together releasing the locking mechanism and gently withdraw the vessel dilator and "J" guidewire, leaving the sheath in place. (Refer to figure 15)
- Warning:** Hold thumb over exposed orifice of sheath to prevent air aspiration. The risk of air aspiration is reduced by performing this part of the procedure with the patient performing the Valsalva maneuver. (Refer to figure 16)
- Insert catheter into lumen of sheath and advance to desired position in vessel. (Refer to figure 17)
- Verify catheter tip location radiographically. **Warning:** This is not a right atrium catheter. Avoid positioning the catheter tip in the right atrium. Placement or migration of the catheter tip into the right atrium may cause cardiac arrhythmia, myocardial erosion or cardiac tamponade. The risk of these potential complications may be more likely in neonatal patients. Preferred location of the catheter tip is at the junction of the superior vena cava and the right atrium.
- Grasp the two handles of the peel-apart sheath and pull outward and upward at the same time. (Refer to figure 18)
- Peel the sheath away from the catheter completely. Make sure the catheter is not dislodged from vessel as sheath is removed.
- Unclamp catheter and withdraw blood through the lumen(s) to insure patency before closing the skin at the venipuncture site. If catheter is not patent, adjust catheter at curvature point to relieve possible restriction. Irrigate catheter lumen(s) with 10ml of sterile normal saline to clear catheter of blood. Instill sterile heparinized saline per lumen to create heparin lock. Clamp catheter.
- Attach injection cap(s) or connect to intravenous fluid source.
- Close the skin at the venipuncture site as necessary, taking care not to damage the catheter.
- Suture catheter at exit site. (Avoid nicking catheter with suture needle.)
- Secure catheter at exit site with a sterile dressing. The external segment of the catheter should be coiled and taped. Avoid tension on the external segment to prevent dislodging the catheter.

Catheter Removal

After tissue grows into the **SureCuff*** Tissue Ingrowth Cuff (2 to 3 weeks), catheters can be removed from the subcutaneous tunnel using one of several methods. The method used will depend upon physician preference and the amount of tissue/cuff ingrowth that is present. The catheter can usually be removed by traction on the external segment (see #1 below) if it is not sutured internally at the cuff or vessel insertion site. Surgical removal (see #2 below) may be necessary to prevent breaking the catheter if the catheter does not dislodge easily with traction or if there is no definite suture site information.

Warning: You should not feel any resistance when withdrawing the catheter from the vein. If you do encounter resistance, this may indicate that the catheter is being pinched between the clavicle and first rib (the "pinch-off" sign). Do not continue pulling against resistance as this may cause catheter breakage and embolism. Free up the resistance (e.g. by repositioning the patient) before proceeding further.

- Traction Removal**
Pull the catheter external segment downward in a straight line away from the exit site with a series of gentle tugs. When separation of the cuff from the surrounding tissue and/or catheter occurs, there will be a "break-away" feeling. Continue to pull gently on the catheter to complete the removal. Apply pressure to the catheter/vein insertion site as needed to control bleeding. If the cuff remains in the subcutaneous tissue, dissect it out through a small incision utilizing local anesthesia.
- Surgical Removal (using aseptic technique)**
 - Locate the position of the cuff either by palpation or by observing the position of "dimpling" when traction is applied to the catheter's external segment.
 - Make a short transverse incision at or below the external side of the cuff taking care not to transect the catheter. Reach under the catheter with a curved, smooth-jawed clamp and pull up on the catheter to remove the catheter tip from the vein. **Caution:** Do not grasp the catheter with any instrument that might sever or damage the catheter.
 - Dissect out the cuff. Transect the catheter on the exterior side of the cuff and remove the interior portion of the catheter and cuff through the incision. **Caution:** Do not cut the catheter before removal from vein to avoid catheter embolism.
 - Remove the exterior segment of the catheter by pulling it from the skin exit site.
 - Apply pressure to the catheter/vein insertion site as needed to control bleeding.
 - Close the incision with a suture as needed. Apply antibiotic ointment to incision and skin exit sites and an occlusive dressing to prevent air embolism through the tract.

Catheter care and maintenance procedures are included in the **Hickman*, Leonard* and Broviac*** CV Catheter Nursing Procedure Manual available through **Bard Access Systems** Customer Service, 1-800-545-0890. For outside the U.S., contact your local sales representative or distributor.

References

1. Aitken, D.R. and Minton, J.P. "The Pinch-Off Sign: A Subclavian Catheters", American Journal of Surgery, Vol. 148, Nov. 1984, pp. 633-636.
2. Rubenstein, R.B., Alberty, R.E., et al. "Hickman* Catheter Separation", JPEN, Vol. 9, No. 6, Nov./Dec. 1985, pp. 754-757.
3. Hinke, D.H.; Zandt-Stastny, D.A.; Goodman, L.R.; et al. Pinch-off syndrome: A complication of implantable subclavian venous access devices. Radiology 177: 353-356, 1990.
4. Ingle, Rebecca; Nace, Corinne, Venous Access Devices: Catheter Pinch-off and Fracture, 1993, **Bard Access Systems**

Patient Information - Catheter Care and Maintenance

Catheter Damage

If the catheter or connection is damaged or dislodged during or after surgery, immediately clamp the catheter with an atraumatic catheter clamp or kink and tape it. The catheter should be repaired as soon as possible using the designated **Hickman***, **Leonard*** and **Broviac*** repair kit for that particular catheter size. Instructions are enclosed in the repair kit package and are also available in the **Hickman***, **Leonard*** and **Broviac*** Catheter Nursing Procedure Manual.

Site Care

Supplies you will need:

- Sterile gloves (if required)
- 3 Alcohol swabsticks
- Hydrogen peroxide
- Sterile cotton-tipped applicators
- 3 Povidone iodine swabsticks
- Povidone iodine ointment packet
- Tape
- 1 Sterile cover dressing (transparent or tape)
- 1 Alcohol wipe
- 1 Sterile 2 in. x 2 in. (5 cm x 5 cm) gauze dressings
- 1 Sterile pre-cut 2 in. x 2 in. (5 cm x 5 cm) gauze dressings

1. Clean the work surface by wiping with a paper towel that has been moistened with alcohol. Wipe dry or allow to air dry. Then place supplies on the cleaned surface.
2. Wash your hands thoroughly using warm soapy water. Rinse completely and dry using a clean towel or fresh paper towels.
3. Carefully open the dressing kit, or unwrap supplies, without touching the inside surfaces of the kits or wrappers.
4. Carefully remove the old dressing, starting from the top of the dressing and working downward. Remove the tape or dressing carefully to avoid irritating your skin or pulling on the catheter.


Caution: Do not use scissors or any sharp-edged instruments as they could damage the catheter.

5. Wash your hands again.
6. Do a careful observation of the exit site and the skin around it. If you notice anything unusual, finish the dressing procedure and then call your doctor.
7. If you are instructed to use gloves, put on the pair of sterile gloves following the procedure you were taught.
Be careful to not touch anything except the supplies being used for site care.
8. Carefully clean the catheter exit site with an alcohol swabstick or sterile cotton-tipped applicator, soaked in hydrogen peroxide, starting at the exit site and spiraling outward until a circle at least 8 cm in diameter, has been cleaned. Do not return to the catheter exit site with a swabstick that has touched any skin away from the exit site. (Refer to figure 19)
9. Repeat this step twice using the other two swabsticks. Look at the color of the swabsticks after you have used them for signs of drainage.
10. Repeat step 8 using three of the povidone iodine swabsticks to clean the same skin area again as well as the part of the catheter that will be lying on the cleaned skin.
11. Gently clean the outside of the catheter with the inside surface of an alcohol wipe, starting from the exit site to the catheter connector. You may hold the catheter at the exit site with another alcohol wipe to prevent pulling on the catheter. **Do Not Pull On The Catheter.** (Refer to figure 20)
12. Allow the povidone iodine on the skin to air dry at least two minutes.
13. Apply a small amount of povidone iodine ointment to the exit site (optional). (Refer to figure 21)
14. Place the pre-cut gauze dressing over the ointment at the exit site, fitting it snugly around the catheter. Place the 2 in. x 2 in. (5 cm x 5 cm) gauze over the pre-cut gauze and catheter.
15. Apply the cover dressing (tape or transparent dressing) following the directions in the package as well as instructions from your doctor or nurse.
16. Coil the catheter, check to see that it is not kinked or pinched, and secure it to the chest or dressing with tape. This will prevent pulling of the catheter at the exit site and decrease irritation.
17. Always secure the catheter in such a way that you can easily see the cap end. Your doctor or nurse will help you select the best method to secure the catheter. The type of clothing and normal activity will need to be considered in this selection. You should periodically look at the capped end to be sure it is intact. (Refer to figure 22)
18. During all dressing changes, assess the external length of the catheter to determine if migration of the catheter has occurred. Periodically confirm catheter placement, tip location, patency and security of dressing.

Clamping the Catheter

Selection of the catheter clamp is very important since the catheter is vital to your care. The wrong clamp can damage the catheter. Follow these three rules for clamping:

1. Use only smooth-edged clamps.
2. Always clamp the catheter over the reinforced clamping sleeve or tape tab, as instructed by your nurse. Never clamp over the reinforced segment directly adjacent to the connector. (see diagram)
3. Follow the directions of your doctor or nurse regarding when to clamp.

Most Hickman* and Broviac* catheters come with pre-attached clamps and reinforced clamping sleeves.

When should you clamp?

Your doctor or nurse may instruct you to clamp your catheter whenever it is not being used. The catheter is filled with heparin and is capped and will protect you from having any problems, but the clamp can be another safety measure. You should always clamp your catheter whenever it is opened to the air, such as during catheter cap changes, or when connecting intravenous infusions to your catheter. **Always have a spare clamp available.**

Flushing the Catheter and “Heparin Lock” Procedure

Supplies you will need:

- Alcohol or povidone iodine wipe.
- 10 ml syringe with attached 1 inch (2.5 cm) needle filled with 2.5 ml of heparin, prepared for use
- Clamp
- Tape

The steps in the procedure are:

1. Collect your supplies in a convenient place.
2. Wash your hands thoroughly.
3. Remove the tape that is around the injection cap.
4. Clean the cap with an alcohol or povidone iodine wipe. If you use the iodine wipe, allow the cap to air dry for two minutes -- be sure not to touch the cap during this time. Do not blow on the area or allow the clean cap to dangle since this increases the chance of contamination of the area with germs.
5. Remove the needle cover and carefully insert the needle into the center of the catheter injection cap. (Refer to figure 23)
6. Release the clamp.
7. Inject the heparin into the catheter. As you inject the last 0.5 ml of heparin solution, withdraw the needle from the injection cap. **If you are flushing the catheter of a child, do not flush too rapidly because the child's circulatory system is small and sensitive to rapid changes in volume and pressure.**
8. Remove the needle from the injection cap. Discard the syringe and needle in a biohazard container.
9. Retape the cap as outlined in the injection cap change procedure.

If you have a multi-lumen catheter, use a separate syringe to flush each lumen with sterile heparin solution. **Your doctor or nurse will give you additional information for the care of multi-lumen catheters.**

Changing the Injection Cap

Supplies you will need:

- Sterile injection cap.
- Alcohol or povidone iodine wipe.
- Catheter clamp
- Tape

The procedure to change the cap:

1. Wash your hands thoroughly.
2. **Be sure the catheter is securely clamped over the reinforced sleeve or tape tab.**
3. Open the package of the new injection cap and prepare according to your instructions. Be sure the cap does not touch the outer surface of the package.
NOTE: You may need to pre-fill the injection cap with heparin if it is a long cap with significant air space. Your doctor or nurse will teach you this additional procedure.
4. Remove the old tape from around the cap by unpeeling the tape. **NEVER** attempt to cut the tape with scissors as you may damage the catheter.
5. Using an alcohol or povidone iodine wipe, clean around the place where the cap is connected to the catheter. Allow to air dry. (Refer to figure 24)
6. **While holding the catheter connector below the level of your heart,** unscrew the old cap and discard. (The fluid level in the catheter will drop part-way into the catheter if the connector is held above the level of your heart.)
7. Pick up the new cap only by the top and remove the sterile tip protector. Attach the new cap by firmly screwing it onto the catheter connector.
8. Cut a 5 cm piece of tape and make tabs on each end by folding back 1 cm. Apply the sticky part of the tape around the connection of the cap and catheter and fasten securely.
Press ends of the tape together. The tabs on the end of the tape will enable you to remove it very easily.
9. Follow the directions your doctor or nurse has given you regarding whether to leave the clamp in place.

Bard Access Systems, Inc. warrants to the original purchaser that this product will be free from defects in material and workmanship for a period of one (1) year from the date of purchase. If this product proves to be so defective, purchaser may return same to Bard Access Systems, Inc. for repair or replacement, at Bard Access Systems, Inc.'s option. All returns must be authorized in advance in accordance with Bard Access Systems, Inc.'s Returned Goods Policy found in its then current Price List. The liability of Bard Access Systems, Inc. under this limited product warranty does not extend to any abuse or misuse of this product or its repair by anyone other than an authorized Bard Access Systems, Inc. representative.

THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED (INCLUDING, WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THE LIABILITY AND REMEDY STATED IN THIS LIMITED PRODUCT WARRANTY WILL BE THE SOLE LIABILITY OF BARD ACCESS SYSTEMS, INC. AND REMEDY AVAILABLE TO PURCHASER FOR THIS PRODUCT, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE) OR OTHERWISE, AND BARD ACCESS SYSTEMS, INC. WILL NOT BE LIABLE TO PURCHASERS FOR ANY INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF ITS HANDLING OR USE.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issued or revision date for these instructions is included for the user's information. In the event two years have elapsed between this date and product use, the user should contact **Bard Access Systems** to see if additional product information is available.

Revised Date: April 2011

***Bard, Broviac, Hickman, Intro-Eze, Leonard** and **SureCuff** are trademarks and/or registered trademarks of C. R. Bard, Inc. All other trademarks are the property of their respective owners.

© 2011 **C. R. Bard, Inc.** All rights reserved.

BARD

Bard Access Systems, Inc.

605 North 5600 West
Salt Lake City, Utah 84116 USA
801-522-5000

Customer Service: 800-545-0890
Clinical Information: 800-443-3385
www.bardaccess.com

Bard Limited

Forest House, Brighton Road
Crawley, West Sussex
RH11 9BP UK

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