ANCORIS

POP REPAIR SYSTEM

Instructions for use
ENGLISH

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INSTRUCTIONS FOR USE

DESCRIPTION

The ANCORIS POP Repair System kit contains:

- 3 TAS anchors (Tissue Anchoring System, Ref: TAS), with the corresponding sutures, manufactured using polypropylene
- 1 Retractable insertion guide (Ref: DPN-MNL), disposable, designed for the placement of the anchors, and provided with a tube in order to protect the tissues during the surgical procedure.
- 1 Knot pusher (Ref: KP), disposable, designed for use during the surgical procedure if necessary, and provided with a protective tube
- 1 Half circle surgical eye needle (Ref: ESN)

All ANCORIS POP Repair System components provided are sterile and ready to use.

INDICATIONS

ANCORIS POP Repair System is indicated for surgical treatment of apical prolapse.

CONTRAINDICATIONS

ANCORIS POP Repair System must not be used in patients:

- Undergoing anticoagulation therapy
- With infectious processes, especially of the genital system or urinary tract
- With vaginal, cervical, or uterine cancer
- During pregnancy
- With a known sensitivity or allergy to polypropylene products
- With pre-existing conditions that pose an unacceptable surgical risk

WARNINGS

ANCORIS POP Repair System must only be used by surgeons and/or operating room staff who are familiarized with the surgical procedures for repairing vaginal defects by the transvaginal approach. Promedon recommends to review the information herein and access an overview of the implant in the official surgical procedure videos (uploaded to the official website) prior to performing the surgical procedure.

Careful patient selection and complete diagnostic evaluation are essential prior to surgery.

Maximum precautions must be taken to avoid contamination of the implant during the surgical procedure. As with all foreign materials, the anchors could exacerbate a preexisting infection.

The surgical procedure must be carried out carefully, avoiding damage to large blood vessels, nerves, and organs. Risks can be minimized by closely observing the local anatomy and by properly inserting the retractable insertion guide.

The described surgical technique is considered to be safe, but in the event of very severe deviations from this and/or extraordinary anatomical variations, perforation or injury may occur to blood vessels, organs, or nerves when passing the retractable insertion guide, and repair surgery may then be necessary.

ANCORIS POP Repair System must be used with caution in patients with:

- Diabetes
- Coagulation difficulties
- Obstruction of the upper urinary tract
- Renal insufficiency
- Autoimmune diseases affecting connective tissue

The patient must be warned that future pregnancies could invalidate the surgical effects of the implant.

ANCORIS POP Repair System components provided are sterile. Do not use the product if the package is open or damaged.

ANCORIS POP Repair System components are designed and intended for SINGLE use only. Therefore, DO NOT REUSE or RE-STERILIZE, since this could decrease the performance of the device and increase the likelihood of cross-contamination.

It is advisable that patients avoid heavy lifting, strenuous exercise (cycling, running, etc.) and sexual intercourse at least during the first four weeks after surgery. The doctor determines when it is appropriate to resume her normal activities.

The patient should contact the surgeon immediately in the event of:

- Dysuria
- Vaginal pain
- Fever
- Serous, bloody or purulent secretions
- Hemorrhages or other problems

PRECAUTIONS

The implant must not be handled with pointed, serrated, or sharp objects since any damage, or tearing can cause subsequent complications.

Maximum precautions must be taken to avoid contamination.

Avoid excessive tension on the implant during insertion.

Maximum precautionary measures must be taken when attaching the TAS to the insertion guide and in order to avoid exerting pressure in the wrong direction during TAS insertion. See "Surgical procedure" section below.

Handling and Storage Precautions

ANCORIS POP Repair System is provided sterile and pyrogen-free. It comes packaged in a cardboard box containing a blister pack and four pouches. The blister pack holds the Retractable Insertion Guide and the Knot pusher. One pouch contains the needle, and the remaining three pouches contain one TAS each and its corresponding suture. **IF ANY OF THE POUCHES AND/OR THE BLISTER PACK ARE DAMAGED, DO NOT IMPLANT.**

Operating room conditions must meet hospital, administrative, or local government requirements.

After use, discard the product and packaging according to hospital, administrative, or local government requirements.

STORAGE

ANCORIS POP Repair System should be stored under the following conditions:

TEMPERATURE: Room temperature

DO NOT USE AFTER THE EXPIRATION DATE SPECIFIED ON THE PACKAGING.

MAGNETIC RESONANCE (MR) ENVIRONMENT

The implant does not affect and is not affected by magnetic resonance (MR) environments.

POSSIBLE COMPLICATIONS

The patient and/or her representatives should be fully informed of the procedure and give consent. The surgeon should discuss with the patient and/or her representatives the possible complications associated with the surgery, anesthesia and device prior surgery.

The implantation of an implant may cause complications related to each patient's specific degree of intolerance to any foreign material implanted in the body. Some complications may require removal of the implant.

Some patients may experience vaginal pain during the initial post-operative period. Treatment with ANALGESICS and ANTI-INFLAMMATORY DRUGS may be sufficient to relieve pain.

Other complications reported with this or other similar implants include:

- Infection
- Adhesion formation
- Vaginal pain, discomfort, irritation
- Purulent, serous or bloody discharge
- Inflammation
- Injuries to blood vessels or nerves
- Presence of vaginal fistula
- Bladder instability
- Bowel problems
- Recurrence of prolapse

Post-operative formation of fibrous tissue around the implant is a normal physiological response to the implantation of a foreign material.

In the event of vaginal exposure of the implant, generally due to infection, it may be necessary to partially remove it.

PROMEDON requires surgeons to report all complications associated with the use of ANCORIS POP Repair System to the company or to the distributor.

PATIENT INFORMATION

The surgeon is responsible for informing the patient and/or her representatives before surgery about the possible complications related to the implantation of the ANCORIS POP Repair System.

The patient should be warned that future pregnancies could invalidate the surgical effects of implanting ANCORIS POP Repair System. The patient should avoid heavy lifting, strenuous exercise (cycling, running, etc.) and sexual intercourse for at least for 4 weeks after surgery. The doctor determines when it is appropriate to resume her normal activities.

The patient should be advised that ANCORIS POP Repair System is a permanent implant, and any complication associated with the implant may or may not require additional surgery to correct the complication.

The patient should contact the surgeon immediately in the event of:

- Dysuria (pain or difficulty with urination)
- Vaginal pain
- Fever
- Serous, bloody or purulent discharges
- Hemorrhages or other difficulties
- Urinary obstruction
- Bowel problems

SURGICAL PROCEDURE

ANCORIS POP Repair System can be implanted by two different approaches, the anterior (paravesical) and the posterior (pararectal) approach.

The decision of which approach is to be preferred depends on the preconditioning of the patient and type of pelvic floor defect:

- In case of a uterine preservation the anterior (paravesical) approach is recommended.
- In case of absence of uterus the approach depends on the type of concurrent defect.
 - In case that an anterior repair needs to be performed in the same surgery it is recommended to implant ANCORIS POP Repair System also via the anterior (paravesical) approach.
 - In case that a posterior repair needs to be performed in the same surgery it is recommended to implant ANCORIS POP Repair System also via the posterior (pararectal) approach.

Anterior Approach:

The description of the technique is summarized in the following steps:

- a. Patient should be in dorsal lithotomy position with legs raised and bent, under spinal, local or general anesthesia. The administration of prophylactic therapy with antibiotics should be considered, according to hospital procedures.
- b. Insert a size 12 or 14 Foley catheter in the urethra.
- c. Make a lengthwise incision along the anterior vaginal wall, stopping 2 cm before the cervix or vaginal vault.
- d. Perform a lateral blunt dissection toward the ischial spine, then identify the coccygeal muscle and the sacrospinous ligament on the right side. The same procedure is performed on the left side.
- e. Identify the right ischial spine as a point of reference and determine the position, size and thickness of the sacrospinous ligament.
- f. To prime the ANCORIS POP Repair System, the tip of the RIG needs to be moved out by pressing the trigger on the handle down and shoving it forward.
- g. Once the tip of the RIG is fully moved out, attach the anchor to the tip by firmly sliding it on to the tip. The TAS is fully connected when the metal tip of the RIG is clearly visible on the top of the TAS. Please note during the whole maneuver that the spikes of the TAS are sharp and may damage gloves.
- h. Optionally shove the protective tube over the RIG and make sure that the sutures of the TAS are running inside the protective tube.
- i. Carefully guide the TAS with the RIG downwards to the sacrospinous ligament by guidance of the index finger.

j. Insert the TAS to the top side of the sacrospinous ligament 2 to 2.5 cm medial from the ischial spine. Figure 1 shows the correct direction to apply pressure when inserting the TAS. The surgeon should use the index finger to touch and identify the ligament and to guide the retractable insertion guide to its correct implant location. The TAS should be placed bilaterally, one in each sacrospinous ligament.


Figure 1

k. Once the TAS have been correctly placed, two stitches are placed at the paracervical ring or vaginal vault with the sutures coming from the TAS and the half circle surgical eye needle provided. The depth of suturing in the attachment points should be deep enough to avoid later tearing.

l. The incision made in the anterior vaginal wall is closed. At this point, the vaginal apex is guided down bilaterally with the help of the index finger and 3 running knots and 2 blocking knots on each side towards the sacrospinous ligaments, avoiding excessive tension.

m. Standard closure of the vaginal wall is completed.

n. Final antisepsis is performed. Digital rectal examination and placement of a vaginal tampon is recommended.

Note: The plastic knot pusher and its protective tube are included in the kit in order to facilitate the knotting procedure of the TAS if the surgeon considers it advisable.

Posterior Approach:

The description of the technique is summarized in the following steps:

a. Patient should be in dorsal lithotomy position with legs raised and bent, under spinal, local or general anesthesia. The administration of prophylactic therapy with antibiotics should be considered, according to hospital procedures.

b. Insert a size 12 or 14 Foley catheter in the urethra.

c. Make a lengthwise incision along the posterior vaginal wall, beginning 2 cm below the vaginal apex.

d. Perform a lateral blunt dissection towards the ischial spine, then identify the ischial spine, the coccygeal muscle and the sacrospinous ligament on the right side. The same procedure is performed on the left side.

e. Identify the right ischial spine as a point of reference and determine the position, size and thickness of the sacrospinous ligament.

f. To prime the ANCORIS POP Repair System, the tip of the RIG needs to be moved out by pressing the trigger on the handle down and shoving it forward.

g. Once the tip of the RIG is fully moved out, attach the anchor to the tip by firmly sliding it on to the tip. The TAS is fully connected when the metal tip of the RIG is clearly visible on the top of the TAS. Please note during the whole maneuver that the spikes of the TAS are sharp and may damage gloves.

h. Optionally shove the protective tube over the RIG and make sure that the sutures of the TAS are running inside the protective tube.

i. Carefully guide the TAS with the RIG downwards to the sacrospinous ligament by guidance of the index finger.

j. Insert the TAS to the top side of the sacrospinous ligament 2 to 2.5 cm medial from the ischial spine. Figure 1 shows the correct direction to apply pressure when inserting the TAS. The surgeon should use the index finger to touch and identify the ligament and to guide the retractable insertion guide to its correct implant location. The TAS should be placed bilaterally, one in each sacrospinous ligament.


Figure 1

- k. Once the TAS have been correctly placed, two stitches are made upon the vaginal apex bilaterally with the sutures coming from the TAS and the half circle surgical eye needle provided. The depth of suturing in the attachment points should be deep enough to avoid later tearing.
- l. The incision made in the posterior vaginal wall is closed halfway. At this point, the vaginal apex is guided down bilaterally with the help of the index finger and 3 running knots and 2 blocking knots towards the sacrospinous ligaments, avoiding excessive tension.
- m. Standard closure of the vaginal wall is completed.
- n. Final antisepsis is performed. Digital rectal examination and placement of a vaginal tampon is recommended.

Note: The plastic knot pusher and its protective tube are included in the kit in order to facilitate the knotting procedure of the TAS if the surgeon considers it advisable.

Postoperative care and therapy are at the surgeon's discretion.

In the event that removal of the implant is required, please note:

If implant removal is necessary due to pain, we recommend trying to cut all the sutures fixating the uterus or vaginal vault first.

If nerves or vessels behind the sacrospinous ligaments are persistently affected after 6-8 weeks, the removal of the TAS from the sacrospinous ligament might be necessary. The risks of this procedure may be higher than the benefits resulting from this removal, so each case should be assessed and decided at the surgeon's discretion.

SYMBOLS USED IN LABELS

CATALOG NUMBER

LOT NUMBER

SERIAL NUMBER

CAUTION

EXPIRY DATE

DO NOT RE-USE

MANUFACTURE DATE

STERILE. STERILIZATION METHOD: ETHYLENE OXIDE

MANUFACTURER

AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

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