



# Composite

**BREAST AUGMENTATION TECHNIQUE**

By Doctor Eric Auclair, plastic surgeon



Sebbin  
PARIS

The partner your body can trust



## The philosophy behind

### COMPOSITE AUGMENTATION

The concept underlying the composite breast augmentation technique is simple: to combine the breast volume augmentation capacities of implants with the remodelling capacities of fat injections. It provides an extremely natural result, even in very thin patients. A reconstructed breast is very similar to a natural breast as the implant is placed in a pre-muscular pocket, in place of the normal mammary gland. Axillary access is also natural because there is no scar on the breast.



Breast implant

+

Lipofilling

for a **natural** result


# The birth

## OF COMPOSITE AUGMENTATION

The first procedure combining breast implants and a subcutaneous fat grafting injection was towards the end of 2006 in a very thin patient with premuscular implants visible around the neckline. At the time, the recommended technique involved replacing implants in a new retromuscular pocket at the expense of a quite serious and painful procedure.

As Dr Delay <sup>[1, 2]</sup> and Dr Coleman <sup>[3]</sup> had just published very promising results for the correction of the after-effects of breast preservation treatment by breast lipomodelling and remodelling, we therefore took the decision to apply the same principles to camouflaging breast implants.

After a few months of use in secondary cases, which enabled us to validate the technique, we decided to apply lipomodelling at the same time as primo-implantation, thereby leading to the birth of the composite breast augmentation technique (CBA) <sup>[4, 5]</sup> which has since been the subject of numerous publications <sup>[6,7,8,10]</sup>.



# The surgical technique

## OF COMPOSITE AUGMENTATION

1

### PATIENT SELECTION

The technique is particularly suited to thin patients, for whom the thickness of the Pinch test around the neckline area is less than 2 cm.

2

### FAT PREPARATION

This starts with liposuction performed using a tumescent technique and a multi-perforated 3mm cannula. Fat is collected in a Fat Washer where it is decanted prior to being retrieved in 10mL syringes for re-injection.

3

### IMPLANTATION

With a pre-axillary access route through an L-shaped incision at the top of the retrofascial, the outer edge of the pectoralis major muscle is reached. Then retrofascial dissection of the implantation pocket is started, guided by the view made possible through the use of specific instruments described later. The purpose of this dissection is to obtain a bloodless, symmetric pocket with the exact dimensions of the implant to be fitted.

# The surgical technique

## OF COMPOSITE AUGMENTATION (CONT'D)

4

### LIPOMODELLING

Once the implant has been fitted, further lipofilling is performed to cover the edges of the implant in areas where it can be noticed, particularly in thin patients. The injection is given using an olive tip cannula to reduce the risk of perforating the implant. It is 15 cm long and 1.2 mm wide. The injection is given from an incision made in the areola, the fat is injected into the subcutaneous space through the Cooper ligaments, close to the subcutaneous vascular network.

5

### POSTOPERATIVE CARE

A sticky bandage is placed on the lipofilling areas at the end of the procedure to facilitate grafting; a medical bra is worn for 3 weeks to maintain the position of the implants.

# Ancillary equipment


## ADAPTED

In order to perform dissection of the prosthetic pocket properly and to control haemostasis, appropriate equipment should be used, and it should be sufficiently long to reach all the corners of the pocket.



This is composed of:

- ✓ Monopolar 35 cm long haemostatic forceps with a curved extremity.
- ✓ L-shaped spreader equipped with cold light lamp and suction. To ensure proper efficacy the blade must be 20 cm long and 3 cm wide.
- ✓ A 15 cm long coagulation tip.



IMPLANT

Which one?






# The choice

## OF IMPLANTS

We only use round implants because adding fat above the upper edge of the implant provides a very natural shape, comparable to that obtained with an anatomical implant, 'anatomization' of the round implant is referred to, without having the risk of rotation which remains one of the drawbacks of anatomical implants.



The round implant is selected based on 2 criteria:

- ✓ The diameter of the implantation base of the breast which must be more or less respected,
- ✓ The projection which will define the size of the cup that is obtained.

## Initial results

Several studies were published from our clinical series: the first publication in 2013<sup>(5)</sup> reported the results of 197 cases, this will be replicated in 2015<sup>(6)</sup> focusing on primary cases. The results of these 2 studies are very encouraging, including capsular contracture less than 1%. In 2014, Prof. Cunningham's team published an article<sup>(9)</sup> on the results of mammograms performed in 52 of our patients which reported that all the examinations are classified BI-RADS 1 or 2 and that none of the images analysed required further examination or biopsy.



**CAPSULAR  
CONTRACTURE**



# For further details

## SEE THE BIBLIOGRAPHY

1. Delay E. Lipomodelling of the reconstructed breast. In: Spear SE, editor. Surgery of the breast: principles and art. 2nd Ed., Philadelphia: Lippincot Williams and Wilkins; 2005.
2. Delay E, Gosset J, Toussoun G, Delaporte T, Delbaere M. Efficacité du lipomodelage pour la correction des séquelles du traitement conservateur du sein [Efficacy of lipomodelling for the correction of after-effects of breast preservation treatment]. *Ann Chir Plast Esthet* 2008;53:153–68.
3. Sydney R Coleman, and A. P. Saboeiro. «Fat Grafting to the Breast Revisited: Safety and Efficacy» *Plast Reconstr Surg* 119, no.3 (2007): 775-785; discussion 786-777.
4. Auclair E., Apport du lipomodelage extraglandulaire dans les augmentations mammaires à visée esthétique (note technique) [Contribution of extra-glandular lipomodelling in breast augmentations for aesthetic reasons]. *Ann Chir Plast Esthet.* (2009); 54; 491-495.
5. Auclair E., Blondeel Ph., Del Vecchio D., Composite Breast Augmentation: Soft tissue planning using implant and fat. *Plast Reconstr Surg*, Sept 2013, Volume 132- Issue 3 - p 558-568.
6. Auclair E, Anavekar N. Combined use of implant and fat grafting for breast augmentation. *Clin Plast Surg.* 2015;42: 307–314, vii.
7. N Kerfant, M.D., AS Henry, M.D, W Hu, M.D., Ph.D., A Marchac, M.D, E Auclair, M.D. Subfascial Primary Breast Augmentation with Fat Grafting: A Review of 156 Cases. *Plast Reconstr Surg*, May 2017, Vol 139 Issue 5, 1080e-1085e.
8. Auclair E., Marchac Alexandre, Book chapter: Composite Breast Augmentation in Atlas of Cosmetic breast surgery published by Pr. Paulo Persichetti, p72-78. 2017, Edizioni Minerva Medical.
9. Jennifer A. Cameron, MD, MPH; Eric Auclair, MD; Michael Nelson, MD; Bruce L. Cunningham, MD, MSc. Radiologic Evaluation of Women Following Cosmetic Breast Augmentation with Silicone Implants and Fat Grafting. *Plast Reconstr Surg*, 2014, Vol. 134 Issue 4S-1, 91-92.
10. Auclair E., Composite breast augmentation for primary and secondary cases, Book chapter in *Fat injection: From filling to regeneration*, edited by S. Coleman, R. Mazzola, L. Pu, 2018. Published by CRC Press.



The partner your body can trust

[www.sebbin.com](http://www.sebbin.com)

GRUPE SEBBIN SAS

39 Parc d'Activités des Quatre Vents - CS80010 Boissy l'Aillierie - 95651 CERGY-PONTOISE CEDEX

Tel.: + 33 1 34 42 13 28 - Fax: + 33 1 34 42 16 88