

SYSTEM FLOW CHART


MAQUET
GETINGE GROUP

VENTILATION
PATIENT CONNECTION



PATIENT CONNECTION VENTILATION

Note! Not all products listed here are available in all countries. For details contact your sales representative.



Line colors:

= Tidal volume (Adult): 100 - 4000 ml

= Tidal volume (Infant): 2 - 100 ml, 10 - 12 mm patient tubes
25 - 350 ml, 15 mm patient tubes



= All tidal volumes (Universal)

= Not included in the patient breathing system



STARTER KIT CONSUMABLES

One selected kit may be ordered with every ventilator

22 mm patient tubes

Starter kit ventilation adult, 22 mm (for passive humidification)
66 96 807

- Including:
 - 1 pc of Breathing system for HME, disposable 22 mm (20 pcs: 66 91 957).
 - 1 pc of Servo Humidifier 162 (100 pcs: 64 19 357).
 - 1 pc of Servo Duo Guard viral/bacterial filter (60 pcs: 66 71 775).
 - 1 pc of NIV reusable face mask, Medium (10 pcs: 66 91 975).
 - 1 pc of NIV face mask, disposable, Size L (20 pcs: 66 91 985).

15 mm patient tubes

Starter kit ventilation pediatric, 15 mm (for passive humidification)
66 96 705

- Including:
 - 1 pc of Breathing system for HME, disposable 15 mm (20 pcs: 66 91 960).
 - 1 pc of Servo Humidifier 161 (100 pcs: 64 19 340).
 - 1 pc of Servo Duo Guard viral/bacterial filter (60 pcs: 66 71 775).
 - 1 pc of NIV reusable face mask, Small (10 pcs: 66 91 972).
 - 1 pc of NIV face mask, disposable, Size S (20 pcs: 66 91 980).

10 mm patient tubes (Not with SERVO-i and SERVO-s)


Starter kit ventilation neonatal, 10 mm (for active humidification)
68 80 511

- Including:
 - 2 pcs of Breathing system for dual heated active humidification 10 mm, insp. limb 2 m, disposable (10 pcs: 66 93 742).
 - 2 pcs of Servo Duo Guard viral/bacterial filter (60 pcs: 66 71 775).

HEAT AND MOISTURE EXCHANGE (HME)


Alt. 1: Passive humidification (HME)

Alt. 2: Active humidification



PATIENT BREATHING SYSTEM

Set of patient breathing system



**Alt. 2: Active humidification
Set of patient breathing system**

Reusable

Disposable

For single heated active humidification (MR 850)

Fisher & Paykel kit of accessories MR 850, reusable
65 95 057

- Including:
 - Heater wire with probe housing.
 - Probe housing adult.
 - Draw wire.
 - Cable clip adult, 3 pcs.
 - Cable clip infant, 3 pcs.
 - Connection cable for reusable tube heater.
 - Tube joint 22 mm out. diam., 15 mm int. diam.

Without heater wire

22 mm patient tubes

Breathing system 22 mm (Heater humidifier, insp. limb 1.6 m, Silicone)
66 95 165

- Including:
 - Patient tube 30 cm (22 mm conn).
 - 4 x Patient tube 60 cm (22 mm conn).
 - Y-piece 22 mm conn.
 - Angled adapter 22/15 mm conn.with luer-lock port.
 - Adapter 22 mm conn.
 - Angled silicone cuff 22 mm conn.
 - 2 x Water trap.

10 mm patient tubes

Breathing system 10 mm (Heater humidifier, insp. limb 1.6 m, Silicone)
66 95 167

- Including:
 - Patient tube 30 cm (10 mm conn).
 - Patient tube 30 cm (22/15 mm conn).
 - 2 x Patient tube 60 cm (22 mm conn).
 - 2 x Patient tube 60 cm (22/10 mm conn).
 - Y-piece 15/10 mm conn.
 - Angled adapter 22/15 mm conn.with luer-lock port.
 - Adapter 22 mm conn.
 - Adapter 15 mm conn.
 - Adapter 10 mm conn.
 - Angled silicone cuff 22 mm conn.
 - 2 x Water trap.

22 mm patient tubes

Breathing system 22 mm for single heated active humidification, insp. limb 1.5 m, Hytrel
66 91 962

- Including:
 - Patient tube 30 cm (22 mm conn).
 - 2 x Patient tube 60 cm (22 mm conn).
 - Patient tube with temp port 110 cm (22 mm conn).
 - Y-piece 22 mm conn.
 - Angled adapter 22/15 mm conn.with luer-lock port.
 - Adapter 22 mm conn.
 - Water trap.

15 mm patient tubes

Breathing system 15 mm for single heated active humidification, insp. limb 1.5 m, Hytrel
66 92 155

- Including:
 - Patient tube 30 cm (22/15 mm conn).
 - Patient tube 60 cm (22/15 mm conn).
 - Patient tube 60 cm (22 mm conn).
 - Patient tube with temp port 110 cm (22/15 mm conn).
 - Y-piece 22/15 mm conn.
 - Angled adapter 22/15 mm conn.with luer-lock port.
 - Adapter 22/15 mm conn.
 - Adapter 15 mm conn.
 - Water trap.

10 mm patient tubes

Breathing system 10 mm for single heated active humidification, insp. limb 1.5 m, Silicone
66 92 162

- Including:
 - Patient tube 30 cm (10 mm conn).
 - Patient tube 30 cm (22/15 mm conn).
 - Patient tube 60 cm (22 mm conn).
 - Patient tube 60 cm (22/10 mm conn).
 - Patient tube with temp port 110 cm (22/10 mm conn).
 - Y-piece 15/10 mm conn.
 - Angled adapter 22/15 mm conn.with luer-lock port.
 - Adapter 15 mm conn.
 - Adapter 10 mm conn.
 - Water trap.


Breathing system 10 mm for single heated active humidification, insp. limb 1.5 m, Hytrel
66 92 157

- Including:
 - Patient tube 30 cm (10 mm conn).
 - Patient tube 30 cm (22/15 mm conn).
 - Patient tube 60 cm (22 mm conn).
 - Patient tube 60 cm (22/10 mm conn).
 - Patient tube with temp port 110 cm (22/10 mm conn).
 - Y-piece 15/10 mm conn.
 - Angled adapter 22/15 mm conn.with luer-lock port.
 - Adapter 15 mm conn.
 - Adapter 10 mm conn.
 - Water trap.

Breathing system 22 mm for single heated active humidification, insp. limb 1.5 m, Silicone
66 91 965

- Including:
 - Patient tube 30 cm (22 mm conn).
 - 2 x Patient tube 60 cm (22 mm conn).
 - Patient tube with temp port 110 cm (22 mm conn).
 - Y-piece 22 mm conn.
 - Angled adapter 22/15 mm conn.with luer-lock port.
 - Adapter 22 mm conn.
 - Water trap.

(Continue to next page)




Alt. 1: Passive humidification (HME)

Alt. 2: Active humidification (F&P MR 850)

PATIENT INTERFACE

Miniflow Patient interface system
(a complete system for providing
NIV NAVA and Nasal CPAP)



Non Invasive Ventilation Masks

Reusable

NIV reusable face mask

Small, 10 pcs/pkg (small headgear size)	66 91 972
Medium, 10 pcs/pkg (standard headgear size)	66 91 975
Large, 10 pcs/pkg (standard headgear size)	66 91 977

• Each piece including:
- 1 NIV reusable face mask.
- 1 headgear.

Disposable

NIV face mask, disposable

Size S, 20 pcs/pkg	66 91 980
Size M, 20 pcs/pkg	66 91 982
Size L, 20 pcs/pkg	66 91 985

• Each piece including:
- 1 NIV disposable face mask.
- 1 headgear.

(Continue to next page)


(Continued from previous page)

PATIENT BREATHING SYSTEM OPTIONS

Reusable



Tube connector 22-15 mm, angled, 5 pcs/pkg
90 37 730




Flex tube 15 mm Hytrel
66 92 455
• Including:
- Adapter 15 mm conn.

VIRAL/BACTERIAL FILTER



Servo Duo Guard viral/bacterial filter, 60 pcs/pkg
66 71 775

Specially designed for use during nebulization.




Servo Guard viral/bacterial filter, 50 pcs/pkg
65 95 487


Specially designed for use when rainout is present in the expiratory breathing circuit.

EDI CATHETER


With Edi module (Not with SERVO-s)




Edi catheter, 16 Fr, adult, 125 cm, 5 pcs/pkg
66 85 790
• Disposable.
• PHT free.




Edi catheter, 8 Fr, adult, 125 cm, 5 pcs/pkg
66 85 782
• Disposable.
• PHT free.




Edi catheter, 12 Fr, pediatric, 125 cm, 5 pcs/pkg
66 85 787
• Disposable.
• PHT free.




Edi catheter, 8 Fr, pediatric, 100 cm, 5 pcs/pkg
66 85 780
• Disposable.
• PHT free.




Edi catheter, 6 Fr, neonatal, 5 pcs/pkg
49 cm 66 85 775
50 cm 66 85 777
• Disposable.
• PHT free.



Guide wire for Edi catheter, 8 Fr, adult, 125 cm, 5 pcs/pkg
66 71 292



Guide wire for Edi catheter, 8 Fr, pediatric, 100 cm, 5 pcs/pkg
66 71 272




Guide wire for Edi catheter, 6 Fr, neonatal, 49 + 50 cm, 5 pcs/pkg
66 71 270

Y SENSOR

With Y sensor measuring module (SERVO-i only)

With Y sensor module (Not with SERVO-i and SERVO-s)




Y sensor adult, 10 pcs/pkg
66 70 660
• Disposable.




Y sensor neonatal, 10 pcs/pkg
66 70 662
• Disposable.



Y-piece 10 mm angled, sampling port, 5 pcs/pkg
68 80 253
• Reusable.



CO₂/Y sensor adapter, 10 pcs/pkg
68 80 104
• Disposable.



Y sensor pressure line, 1.8 m, disposable, 10 pcs/pkg
66 93 672

Reusable

Disposable



Y sensor neonatal, reusable
66 93 670



Y sensor neonatal, disposable, 5 pcs/pkg
66 93 667

CO₂ MONITORING


With CO₂ analyzer module (Not with SERVO-s)

Reusable


Disposable

Reusable


Disposable




CO₂ airway adapter, adult (reusable)
47 21 796



CO₂ airway adapter, adult (disposable), 10 pcs/pkg
66 84 612



CO₂ airway adapter, neonatal (reusable)
47 21 788




CO₂ airway adapter, neonatal (disposable), 10 pcs/pkg
66 84 610

NEBULIZATION

With Aereoneb

With Servo Ultra Nebulizer

Note! The Servo Ultra Nebulizer has been discontinued. The information in this document applies to devices sold with versions 6.1 and earlier of the SERVO-i Ventilator System



Contents

- Starter kit consumables
- Heat and Moisture Exchange (HME)
- Patient breathing system
- Patient interface
- Patient breathing system options
- Viral/bacterial filter
- Edi catheter
- Y sensor
- CO₂ monitoring
- Nebulization
- Test lung

For local contact:

Please visit our website
www.maquet.com

MAQUET
GETINGE GROUP

Maquet Critical Care AB
SE-171 54 Solna
SWEDEN
Phone: +46 (0) 8 730 73 00
www.maquet.com

Products may be pending regulatory approvals to be marketed in your country. Contact your MAQUET representative for more information.

GETINGE GROUP is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, GETINGE and MAQUET. ArjoHuntleigh focuses on patient mobility and wound management solutions. GETINGE provides solutions for infection control within healthcare and contamination prevention within life sciences. MAQUET specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.