



ITEM	TITLE	120x70 cm	140x70 cm	160x70 cm	180x70 cm	190x70 cm
7	Sengehest komplet *	40-30196-1	40-30196-1	40-30197-1	40-30197-1	40-30197-1
8	Løsdele til sengehset	40-30198-1	40-30198-1	40-30199-1	40-30199-1	40-30199-1
9	Bordplade mobilio	40-30191-1	40-30192-1	40-30193-1	40-30194-1	40-30195-1

40-30201-1	1	Løftesøjle inkl. beslag	1		6
97001655	1	Håndbetjening	1		5
97001724	1	Motorkabel	1		4
97000427	1	Netkabel			3
97001720	1	Kontrolboks CO61	1		2
40-30200-1	4	Hjul ø100 mm. med brems	1		1
PART NUMBER	QTY	TITLE	L	B	ITEM

*Sengehest komplet indeholder nødvendige løsdele til montage.
 *Sengehest er tilbehør til nogen modeller

1										
2										
3										
4										
5										
Side 2 af 2						ROPOX Better ways to better days Ringstedgade 221, 4700 Næstved Tlf.: 55750500 Fax.: 55750550			Tegning nr.: 40-30406-5	

Tegn.	sj	13-11-2019
DS/ISO 2768 m		
Vægt:	N/A	N/A