


- Axila crutch of beechwood, varnished natural
- PU foam padding
- Rubber tip made of natural rubber



5-fold height adjustable
by metal pin




Inner-Ø 20 mm



Tread Ø 35 mm
Article no. 920gr




in the lower part 6-fold
height adjustable by 2
splints and wing nuts



15K

Measure A
Armpit height
110 cm – 135 cm

Measure B
Grip height
39 cm – 49 cm



15G

Measure A
Armpit height
132 cm – 151 cm

Measure B
Grip height
50 cm – 60 cm



MAX. load bearing capacity
up to 150 kg body weight



WEIGHT
1800 g / Pair (Article 15K)
2000 g / Pair (Article 15G)

