

# Incidin 70

## Disinfecting Wipes

Alcohol wipes for the disinfection of non-invasive medical devices leaving surfaces hygienically clean in seconds.

### AN EFFECTIVE ALCOHOL BLEND

Incidin 70 wipes are a blend of 70% denatured ethanol and 5% IPA for the disinfection of non-invasive medical devices.

Tested against a wide range of micro-organisms, Incidin 70 wipes leave surfaces hygienically clean without leaving any residues.

### STRONG WIPE MATERIAL

Incidin 70 wipes use the latest in synthetic wipe innovation. Incidin 70 wipes are super absorbent and incredibly low linting.

### WIPE SPECIFICATION

**Ingredients:** 70% Denatured Ethanol, 5% IPA, Water

**Expiry date:** 36 months

**Shelf life once opened:** 6 months

**Wipe material:** White Tritex™ 18gsm



### INCIDIN 70 EFFICACY TABLE

EN TEST METHODOLOGY	ORGANISM	CONDITIONS	CONTACT TIME
<b>BACTERICIDAL EFFICACY</b>			
EN 13727	Methicillin - resistant Staphylococcus aureus	Dirty	1 minute
	Staphylococcus aureus		
	Enterococcus hirae		
	Pseudomonas aeruginosa		
	Vancomycin - resistant enterococcus faecalis		
EN 13697	Acinetobacter baumannii	Dirty	1 minute
	Pseudomonas aeruginosa		
EN 14561	Escherichia coli	Dirty	1 minute
	Staphylococcus aureus		
	Pseudomonas aeruginosa		
	Acinetobacter baumannii		
	Escherichia coli		
EN 13624	Enterococcus hirae	Clean	5 minutes
	Klebsiella pneumoniae		
	Candida albicans		
	Aspergillus brasiliensis		
EN 14652	Aspergillus niger	Dirty	1 minute
	Candida albicans		
<b>MYCOBACTERICIDAL EFFICACY</b>			
EN 14653	Mycobacterium terrae	Dirty	1 minute
	Mycobacterium avium		
EN 14348	Mycobacterium terrae	Dirty	1 minute
	Mycobacterium avium		
<b>VIRUCIDAL EFFICACY</b>			
EN14476	Influenza A H1N1	Clean	1 minute

### PRODUCT INFORMATION

PRODUCT DESCRIPTION	PRODUCT CODE	WIPE SIZE	WIPE COUNT	CASE SIZE
Incidin 70 Alcohol Wipes	3077450	195mmx210mm	150 Wipes	10 Tubs per case

ALWAYS READ THE LABEL & PRODUCT INFORMATION BEFORE USE.

#### MARKEDSFØRES AV

Puls a.s.  
Strømsveien 344  
1081 Oslo  
+47 23 32 30 00  
www.puls-norge.no

#### PRODUZERT AV

Ecolab Deutschland GmbH  
Ecolab-Allee 1  
40789 Monheim am Rhein  
+49 (0) 2173 599 0

