

SOL-CARE™ Safety Blood Collection Needle

with or without pre-attached holder



Discover our range of safety 'butterfly' needles for blood collection. Each device features a simple safety mechanism that should be activated immediately after blood withdrawal is complete, thereby reducing the risk of needlestick injury.

- Thin walled needle minimises patient discomfort
- Flexible wings comfortably secure device to patient's skin
- DEHP-free tubing is soft, flexible and kink-resistant
- Transparent design allows for clear visualisation of blood flow

With pre-attached holder




**Simple, intuitive safety system
minimises risk of needlestick injury:**

- Three simple ways to activate safety mechanism – using either one or two handed technique
- Audible click tells you when the needle is safely locked within the safety shield

SOL-CARE™ Safety Blood Collection Needle Range

Product configuration	Product code	Needle outer diameter		Needle length		Tubing length	
		Gauge	Mm	Inches	Mm	Inches	Mm
Luer Adapter without Pre-attached holder	110201030009	25G	0.5	3/4	19	12	305
With Pre-attached holder	110201030010						
Luer Adapter without Pre-attached holder	110201030005	23G	0.6	3/4	19	12	305
With Pre-attached holder	110201030006						
Luer Adapter without Pre-attached holder	110201030001	21G	0.8	3/4	19	12	305
With Pre-attached holder	110201030002						
Luer Adapter without Pre-attached holder	110201030011	25G	0.5	3/4	19	7	178
With Pre-attached holder	110201030012						
Luer Adapter without Pre-attached holder	110201030007	23G	0.6	3/4	19	7	178
With Pre-attached holder	110201030008						
Luer Adapter without Pre-attached holder	110201030003	21G	0.8	3/4	19	7	178
With Pre-attached holder	110201030004						



For more information about SOL-CARE™ Safety Blood Collection Needles, please contact your local SOL-MILLENNIUM® representative.

Sol-Millennium Europe Sp. z o.o.
 Twarda 18, PL-00105 Warsaw, Poland
 Tel: +48 22 202 67 23, Fax: +48 22 202 67 19
 info.eu@sol-m.com