



Power. Saline. One.

BAIRD



**THE POWER
OF PURPLE**

Power
Injectable PICC

Power.

The Efficiency of Power Injection for Contrast Enhanced CT Scans and PICC Therapy.


- ▶ **Maximum Injection Rates**
 - allows injection of contrast media for CECT scans at a maximum rate of **5ml/sec**
- ▶ **Easy Identification**
 - catheters are clearly labeled to differentiate them from all non-power injectable PICCs
- ▶ **CVP Monitoring**
 - all POWERPICC SOLO²® catheters are indicated for central venous pressure monitoring
- ▶ **Reverse Taper Design**
 - gently plugs the insertion site and provides superior kink resistance
- ▶ **Patient Satisfaction**
 - improves patient comfort due to decreased needle sticks

4 Fr. Single-Lumen
POWERPICC SOLO²® Catheter

(Also Available in 5 Fr. Single-Lumen)

5 Fr. Dual-Lumen
POWERPICC SOLO²® Catheter


6 Fr. Triple-Lumen
POWERPICC SOLO²® Catheter



Saline.

The Simplicity of Saline Care and Maintenance Through Innovative Valve Technology.

- ▶ **Simplified Care and Maintenance**
 - recommended weekly maintenance is reduced to a single saline flush
- ▶ **Continuum of Care**
 - same flushing protocol for hospital and home health
- ▶ **Fewer Complications**
 - limits risk associated with heparin use (HIT)
- ▶ **Reduced Blood Reflux**
 - compared to open-ended catheters
- ▶ **Reduced Costs**
 - decreases PICC associated heparin costs
- ▶ **Clamp Free Infusion Therapy**
 - valve eliminates the need for clamps and the maintenance challenges associated with needleless injection caps
- ▶ **Superior Clearance**
 - valve clears after recommended saline flush



One.

The Versatility of One PICC Designed to Meet a Variety of Clinician and Patient Needs.

- ▶ **Broad Vascular Access Needs**
 - designed for multiple clinical conditions across the hospital
- ▶ **Multiple Indications**
 - allows IV fluid and blood infusion, as well as power injection, and central venous pressure monitoring
- ▶ **Many Configurations**
 - 4 and 5 Fr. Single-Lumen, 5 Fr. Dual-Lumen, 6 Fr. Triple-Lumen
- ▶ **Variety of Placement Technologies**
 - Interventional Radiology - allows for over-the-guidewire placement

BLUE COLOUR

Symbolises Saline
Care and Maintenance



POWERPICC SOLO²® catheter is a universal catheter solution designed for IV fluid and blood administration, power injection, and saline care and maintenance.

POWERPICC SOLO²® Catheter



Product Code	Basic Nursing Tray	Guidewire Length	Gravity Flow Rate	Gauge
6194118	4 Fr. SL POWERPICC SOLO ² ® Catheter	50cm	1024	18
6195118	5 Fr. SL POWERPICC SOLO ² ® Catheter	50cm	978	18
6295118	5 Fr. DL POWERPICC SOLO ² ® Catheter	50cm	498/498	18/18
6396118	6 Fr. TL POWERPICC SOLO ² ® Catheter	50cm	R560 G308 W308	R17 G19 W19

Basic Nursing Tray Kit Content

No. 11 Scalpel - IV Catheter 20Ga. - Safety Introducer Needle, 21Ga. - Microintroducer 7cm - Paper Measuring Tape Dead end caps (1 per lumen) - STATLOCK® Catheter Stabilisation Device - Syringe 12ml

Product Code	IR Tray	Guidewire Length	Gravity Flow Rate	Gauge
6194335	4 Fr. SL POWERPICC SOLO ² ® Catheter	135cm	1024	18
6194355	4 Fr. SL POWERPICC SOLO ² ® Catheter	70cm		
6195335	5 Fr. SL POWERPICC SOLO ² ® Catheter	135cm	978	18
6195355	5 Fr. SL POWERPICC SOLO ² ® Catheter	70cm		
6295335	5 Fr. DL POWERPICC SOLO ² ® Catheter	135cm	498/498	18/18
6295355	5 Fr. DL POWERPICC SOLO ² ® Catheter	70cm		
6396335	6 Fr. TL POWERPICC SOLO ² ® Catheter	135cm	R560 G308 W308	R17 G19 W19
6396355	6 Fr. TL POWERPICC SOLO ² ® Catheter	70cm		

IR Tray Kit Content

No. 11 Scalpel - IV Catheter 20Ga. - Introducer Needle, 21Ga. - Microintroducer 10cm - Paper Measuring Tape Dead end caps (1 per lumen) - STATLOCK® Catheter Stabilisation Device - Syringe 12ml



- Kink resistant guidewire
- Bendable proximal end
- Length of guidewire offers greater placement versatility



Learn More

For more information on POWERPICC SOLO²®, please contact your Bard Access Systems representative. Or visit us on the web at www.bardaccess.com or www.powerpiccsolo.com

Bard France S.A.S. (France, North Africa) FRANCE
Tel: + 33 1 39 30 58 58 Fax: + 33 1 39 30 58 59
service.commande@crbard.com www.bardfrance.fr

Bard de España S.A. SPAIN
Tel: + 34 93 253 7800 Fax: + 34 93 253 7834
spain.customerservice@crbard.com

C. R. Bard Portugal L.D.A. PORTUGAL
Tel: + 351 21 319 0330 Fax: + 351 21 319 0339
Bard.Portugal@crbard.com

Bard S.p.A. (Italy, Eastern Mediterranean, Middle East) ITALY
Tel: + 39 06 524 93 249 Fax: + 39 06 529 58 52
customerservice.italia@crbard.com www.barditalia.it

Bard Hellas SA GREECE
Tel: + 30 210 96 90 780 1 Fax: + 30 210 96 28 810
Bard.Hellas@crbard.com

Bard Medica S.A. SWITZERLAND
Tel: + 41 44 722 5360 Fax: + 41 44 722 5370
info-bce.access@crbard.com

Bard Czech Republic S.R.O. CZECH REPUBLIC
Tel: + 420 242 408 630 Fax: + 420 242 410 185
info-bce.access@crbard.com

C. R. Bard GmbH GERMANY
Tel: + 49 721 9445 124 Fax: + 49 721 9445 100
info-bce.access@crbard.com

Bard Poland Sp. z o.o. POLAND
Tel: + 48 22 321 0930 Toll Free: 0800 494900
Fax: + 48 22 321 0938 Toll Free: 0800 494901
info-bce.access@crbard.com

Bard Medica S.A. AUSTRIA
Tel: + 43 1 49 49 130 Fax: + 43 1 49 49 130 30
info-bce.access@crbard.com

Bard Limited UNITED KINGDOM
Tel: + 44 1293 527888 Fax: + 44 1293 552428
customers.services@crbard.com www.barduk.com

IRELAND
Tel: + 353 1 2130706 Fax: + 353 1 2130604
customers.services@crbard.com

Bard Norden AB SWEDEN
Tel: + 46 42 38 60 00 Fax: + 46 42 38 60 10
support.nordic@crbard.com

Bard Norway AS NORWAY
Tel: + 47 815 00 266 Fax: + 47 815 00 299
support.nordic@crbard.com

Bard Denmark AB DENMARK
Tel: + 45 49 26 20 08 Fax: + 45 49 26 20 09
support.nordic@crbard.com

Bard Finland Oy FINLAND
Tel: + 358 9 347 30 53 Fax: + 358 9 347 10 07
support.nordic@crbard.com

Bard Medical SA (Pty) Ltd. SOUTH AFRICA
Tel: + 27 11 524 9900 Fax: + 27 86 537 7250
www.bard.co.za

Bard Benelux N.V. BELGIUM
Tel: + 32 14 28 69 52 Fax: + 32 14 28 69 66
custservicebelgium@crbard.com

THE NETHERLANDS
Tel: + 31 30 60 00 572 Fax: + 31 30 60 22 345
custservicedenland@crbard.com

