




VALVE/CHANNEL CLEANING BRUSHES



Handy cleaning of working channels and valves

- ✓ High cleaning efficiency
- ✓ Rounded brush tips
- ✓ Maximum brush density



For flexible Endoscopes



Valve/channel brushes

Handy cleaning of working channels and valves



Features

- ☑ Rounded brush tips for gentle cleaning
- ☑ Brushes made of soft, durable nylon
- ☑ Rounded handle

TWO-SIDED VALVE/CHANNEL CLEANING BRUSH



Colour	Length	Valve brush 1		Channel brush 2		Pieces per box	Manufacturer Product Code
		Ø	Length	Ø	Length		
Blue	15 cm	11 mm	35 mm	5 mm	20 mm	200	DVS-30-35-10-20-05-15

VALVE CLEANING BRUSH WITH HANDLE



Colour	Length	Valve brush		Pieces per box	Manufacturer Product Code
		Ø	Length		
Green	9.5 cm	11 mm	35 mm	200	SVS-30-35-11

CHANNEL CLEANING BRUSH FOR ALBARRÀN LEVER ON DUODENOSCOPES



Colour	Length	Channel brush		Pieces per box	Manufacturer Product Code
		Ø	Length		
Red	8 cm	2 mm	15 mm	200	SES-30-2-15