





VALVE SETS FOR ENDOSCOPE



## VALVE SETS



### Valves and adapters in a practical disposable set

- ✓ Olympus and Fujifilm
- ✓ High safety
- ✓ Constant performance




### Valve sets for endoscope

In a practical disposable set



#### OLYMPUS



Set contents

#### Features

Biopsy valve, air/water valve, suction valve

- ☐ Suitable for Olympus OES 40, EVIS 140, 145, 160, 165, 180, 185, 190, 240, 260, 290
- ☐ Minimises cross-contamination risks
- ☐ Supports constant performance

#### Jet channel adapter

- ☐ Used together with the rinsing hose in connection with a rinsing pump for flushing with rinsing fluids
- ☐ With integrated non-return valve
- ☐ Individually sterile packed

	Pieces per box	Manufacturer Product Code
Biopsy valve, air/water valve, suction valve & jet channel adapter	25	AWV-KIT-01-OLY
Biopsy valve, air/water valve, suction valve	25	AWV-KIT-02-OLY
Air/water valve, suction valve	25	AWV-KIT-03-OLY
Biopsy valve, air/water valve, suction valve, air/water cleaning valve	25	AWV-KIT-04-OLY
Biopsy valve, air/water valve, suction valve, air/water cleaning valve & jet channel adapter	25	AWV-KIT-05-OLY
Air/water valve, suction valve, air/water cleaning valve & jet channel adapter	25	AWV-KIT-06-OLY
Air/water cleaning valve	100	AWV-KIT-07-OLY
Jet channel adapter	100	ADA-JC-OUS

#### FUJIFILM



Set contents

#### Features

- ☐ Suitable for Fujifilm series
- ☐ Minimises cross-contamination risks
- ☐ Supports constant performance

	Stück je Box	Manufacturer Product Code
<b>SERIES 700</b>		
Biopsy valve, air/water valve, suction valve	25	AWV-KIT-01-FJ
Air/water valve, suction valve	25	AWV-KIT-02-FJ
Air/water cleaning valve	100	AWV-KIT-07-FJ
Jet channel adapter	100	ADA-JC-FJ
<b>SERIES 500/600</b>		
Biopsy valve, air/water valve, suction valve	25	AWV-KIT-03-FJ
Air/water valve, suction valve	25	AWV-KIT-04-FJ



Company GmbH



**UpViser Aps**  
Adelgade 45, 1  
8660 Skanderborg, Denmark  
Tel.: +45 6171 7625  
E-Mail: salg@upviser.com  
www.upviser.dk



**Exclusive Partner Europe**  
Body Products GmbH

**mdd Company GmbH**  
Rengglochstr. 38  
CH-6012 Obernau, Switzerland  
Tel: +39 348 0026970  
E-Mail: gb@mddcompany-gmbh.com  
www.mddcompany-gmbh.com