



BITE BLOCKS



Protection and hygiene for patients and endoscopes

- ✓ Ergonomically shaped
- ✓ Soft edges
- ✓ With and without holding straps

**Bite blocks**

Protection and hygiene for patients and endoscopes

**Features**

- Soft edges increase safety and patient protection
- Optional holding straps increase patient comfort
- Single-use articles reduce the risk of cross contamination

ADULT – OXYGEN**Features**

- Comfortable delivery of additional oxygen through integrated connections

Colour	Product description	Type of strap	Opening dimensions	Pieces per box	Manufacturer Product Code
Light blue	No strap	-	25 x 18 mm	100	BB-05-30-OXY
Light blue	With white textile strap	Textile	25 x 18 mm	100	BB-05-30-OXY-T
Light blue	With rubber strap	Rubber	25 x 18 mm	100	BB-05-30-OXY-R

ADULT**Features**

- Large lateral openings optimise control during handling and improve access
- Fits facial contours thanks to ergonomic design

Colour	Product description	Type of strap	Opening size	Pieces per box	Manufacturer Product Code
Silver	No strap	-	25 x 22 mm	100	BB-05-30-30
Silver	With blue textile strap	Textile	25 x 22 mm	100	BB-05-30-30-T
Silver	With rubber strap	Rubber	25 x 22 mm	100	BB-05-30-30-R



Optimal safety

**Bite blocks**

Protection and hygiene for patients and endoscopes

**Features**

- Soft edges increase safety and patient protection
- Optional holding straps increase patient comfort
- Single-use articles reduce the risk of cross contamination

ADULT**Features**

- Large lateral openings optimise control during handling and improve access
- Fits facial contours thanks to ergonomic design

Colour	Product description	Type of strap	Opening size	Pieces per box	Manufacturer Product Code
Yellow	No strap	-	27 x 23 mm	100	BB-01-85-50
Yellow	With white textile strap	Textile	27 x 23 mm	100	BB-01-85-50-T
Yellow	With green rubber strap	Rubber	27 x 23 mm	100	BB-01-85-50-R

PAEDIATRIC**Features**

- Suitable for children and people with small mouths
- Large lateral openings optimise control during handling and improve access
- Fits facial contours thanks to ergonomic design

Colour	Product description	Type of strap	Opening size	Pieces per box	Manufacturer Product Code
Orange	No strap	-	19 x 16 mm	100	BB-03-30-33
Orange	With white textile strap	Textile	19 x 16 mm	100	BB-03-30-33-T
Orange	With rubber strap	Rubber	19 x 16 mm	100	BB-03-30-33-R