



## ENDOSCOPE TRANSPORT TRAY



### Safe methods for transporting processed and used endoscopes

- ✓ High stability
- ✓ Light-weight
- ✓ Stapelbar



### Endoscope transport tray & CC-cover

Safe methods for transporting processed and used endoscopes



The cardboard endoscope transport trays are meant for single-use and are used to individually transport both sterilized and used endoscopes safely. The practical transport trays are sagely sealed in a CC-cover and the endoscope is marked either as 'clean' or 'contaminated'.

#### CARDBOARD ENDOSCOPE TRANSPORT TRAY WITH CC-COVER.

Just turn over the cc-cover to mark whether the endoscope is 'clean' or 'contaminated' (smiley)

Reversible cc-cover with 'clean' and 'contaminated' smileys

Endoscope transport tray

Stackable transport tray

Keeps an endoscope safely in place during transport through the texture of its surface

Keeps the transport tray safely sealed through the elastic band around the cc-cover

- Characteristics**
- Endoscope transport tray**
- ☑ Moulded trays which are also suitable for longer, flexible endoscopes
  - ☑ The structured surface of the tray holds the endoscope in place during transport.
  - ☑ Very stable due to the cardboard's thickness of 900 g/m<sup>2</sup>
  - ☑ Low weight
  - ☑ Light grey cardboard
  - ☑ Conveniently stackable
  - ☑ Compostable
  - ☑ Sterilized by gamma radiation
- Reversible CC-cover**
- ☑ The 'clean' and 'contaminated' smileys are on transparent polyethylene film (PE) and can be changed by turning the cover over.
  - ☑ The elastic band around the CC-cover is flexible and
  - ☑ saves space when rolled up

Description	Dimensions	Packaging unit	Manufacturer Product Code
<b>Cardboard transport tray</b>	52 x 42 x 10 cm (L x W x H)	10	SUT-52-42-001
<b>CC-cover made out of PE</b> with 'clean' & 'contaminated' smileys	70 x 58 cm (L x W) (Suitable for 52 x 42 x 10 cm)	4 x 25	SUT-52-42-LIN
<b>Transport tray &amp; cc-cover set</b>	As described above	10	SUT-52-42-SET



## ENDOSCOPE TRANSPORT AND PROTECTIVE COVERS



### Protection & clear labelling of processed and used endoscopes

- ✓ “Clean” and “contaminated” labels
- ✓ Transparent polyethylene film (PE)
- ✓ Self-closing, flexible elastic band



### Transport bag for endoscope transport trays

Protection & clear marking of processed and used endoscopes



The practical covering system is suitable for safely labelling of all common transport trays and for protecting endoscopes. The transport trays are safely covered with a protective bag and are marked as either 'clean' or 'contaminated' with the cc-transport bag.

#### PROTECTIVE BAG & CC-TRANSPORT BAG

Protective bag

CC-transport-bag with "Clean" & "Contaminated" marking (text version) for turning around

Simply turn around the CC-transport-bag from "Clean" to "Contaminated", if the endoscope is contaminated

**Features**

**Protective bag**

- ✔ To protect the endoscope or the transport tray
- ✔ Saves space when rolled up
- ✔ Self-closing, flexible elastic band

**CC-transport-bag for turning around**

- ✔ "Clean" and "Contaminated" marking on transparent polyethylene film (PE)
- ✔ Saves space when rolled up
- ✔ Self-closing, flexible elastic band
- ✔ Not available in some countries

Length x width	Suitable for transport trays to length x width x height	Packaging unit	Manufacturer Product Code
S, 70 x 58 cm	52 x 42 x 10 cm	4 x 25	STB-01-SET
XXL, 95 x 80 cm	64 x 46 x 9 cm	4 x 25	STB-02-SET



Safe and protected transport



### Transport bag for endoscope transport trays

For labelling sterilized and used endoscopes quickly



The C transport bag with a C-foil cover is suitable for labelling contaminated endoscopes quickly. The 'clean' smiley is covered up by simply turning the C-foil cover over and then the endoscope will be marked as 'contaminated' without the transport bag having to be removed.

#### 'CLEAN' TRANSPORT BAG WITH 'CONTAMINATED' FOIL COVER

Slide endoscope transport tray into C-transport-bag with "clean" marking

Flip over C-foil-flap with "Contaminated" marking, if the endoscope is contaminated

**Features**

**C-transport-bag with C-foil-flap**

- ☞ "Clean" marking on transparent C-transport-bag (PE)
- ☞ C-foil-flap with "Contaminated" marking for quick cover of the "Clean" marking
- ☞ Simple application
- ☞ Saves space when rolled up
- ☞ Self-closing, flexible elastic band
- ☞ Not available in some countries

Length x width	Suitable for transport trays to length x width x height	Packaging unit	Manufacturer Product Code
S, 70 x 58 cm	52 x 42 x 10 cm	4 x 25	SEC-01
XXL, 95 x 80 cm	64 x 46 x 9 cm	4 x 25	BEC-01