

HAEMOSTASIS



Company GmbH



INJECTION NEEDLES

HEMOCLIPS



INJECTION NEEDLES




Application of medications and solutions

- ✓ Automatic needle retraction
- ✓ Needle locking
- ✓ Colour coding

mdd

Injection needles

Application of medications and solutions



Features

- ☐ Reliable needle guidance thanks to automatic needle retraction
- ☐ Rapid selection of the appropriate needle diameter thanks to permanent colour code
- ☐ Smooth puncturing of the target tissue
- ☐ Needle locking with extended or retracted needle
- ☐ Integrated ergonomic handle for comfortable use
- ☐ Individually sterile packed

Needle diameter	Needle length	Outer catheter diameter	Working length	Pieces per box	Manufacturer Product Code
FOR COLONOSCOPY					
19 gauge / 1.0 mm	4 mm	2.4 mm	230 cm	10	IN-23-19-4-230
19 gauge / 1.0 mm	5 mm	2.4 mm	230 cm	10	IN-23-19-5-230
19 gauge / 1.0 mm	6 mm	2.4 mm	230 cm	10	IN-23-19-6-230
22 gauge / 0.7 mm	4 mm	2.4 mm	230 cm	10	IN-23-22-4-230
22 gauge / 0.7 mm	5 mm	2.4 mm	230 cm	10	IN-23-22-5-230
22 gauge / 0.7 mm	6 mm	2.4 mm	230 cm	10	IN-23-22-6-230
23 gauge / 0.6 mm	4 mm	2.4 mm	230 cm	10	IN-23-23-4-230
23 gauge / 0.6 mm	5 mm	2.4 mm	230 cm	10	IN-23-23-5-230
23 gauge / 0.6 mm	6 mm	2.4 mm	230 cm	10	IN-23-23-6-230
25 gauge / 0.5 mm	4 mm	2.4 mm	230 cm	10	IN-23-25-4-230
25 gauge / 0.5 mm	5 mm	2.4 mm	230 cm	10	IN-23-25-5-230
25 gauge / 0.5 mm	6 mm	2.4 mm	230 cm	10	IN-23-25-6-230



HEMOCLIPS



Hemoclip for endoscopic marking and hemostasis

- ✓ Secure clip placement
- ✓ 360° rotation
- ✓ Colour coding

mdd

Heliclip

Hemoclip for endoscopic marking and hemostasis



Features

- ☐ 180 degree opening and stepless 360° rotation of the clip facilitate placement in the tissue
- ☐ Safe and easy release of the clip
- ☐ PTFE tube coating ensures optimum sliding properties
- ☐ Atraumatic tip protects the endoscope from damage
- ☐ Ergonomic handle supports a comfortable application
- ☐ Colour-coded for quick orientation in choosing the appropriate clip length
- ☐ Individually sterile packed
- ☐ Available from March*, June**, July 2020***

Clip length	Opening	Outer catheter diameter	Working length	Pieces per box	Manufacturer Product Code
9 mm***	18 mm	1.8 mm	230 cm	10	HC-23-230-9
12 mm**	24 mm	1.8 mm	230 cm	10	HC-23-230-12
16 mm*	32 mm	1.8 mm	230 cm	10	HC-23-230-16



Company GmbH



UpViser Aps
Adelgade 45, 1
8660 Skanderborg, Denmark
Tel.: +45 6171 7625
E-Mail: salg@upviser.com
www.upviser.dk



Exclusive Partner Europe
Body Products GmbH

mdd Company GmbH
Rengglochstr. 38
CH-6012 Obernau, Switzerland
Tel: +39 348 0026970
E-Mail: gb@mddcompany-gmbh.com
www.mddcompany-gmbh.com