

Syringe Rack

- For temporary depositing of:
 - max. six 50 ml syringes and
 - max. eight 10 ml syringes
- Made of stainless steel
- Autoclavable
- Dimension 230 x 83 x 95 mm



Ordering Information

	RE	VE	
Syringe Rack	610300	1 pc	supplied non-sterile, autoclavable, reusable

© human med AG 04/2011

human med AG
Wilhelm-Hennemann-Str. 9
19061 Schwerin, Germany
Tel: +49 (0) 385 / 3 95 70 - 0
Fax: +49 (0) 385 / 3 95 70 - 1
info@humanmed.com

