

STONE BASKETS

StoneWrapper[®]

The high performer of stone extraction baskets...

„StoneWrapper[®]“



Perfect for removing concretions of different sizes from hard-to-reach areas of the renal pelvis system comfortably, safely and quickly.

- ◆ Perfect hybrid of grabber and stone extraction basket
- ◆ Elastic, fully articulated three-armed nitinol basket
- ◆ Easy insertion into the lower calyx and quick removal of all concretions
- ◆ Available in two sizes and with two opening options

Effective: work efficiently, save time

The innovative StoneWrapper solves the major challenges involved with flexible ureterorenoscopy: fast and safe atraumatic removal of marginal stones from the renal pelvic system.

STONE BASKETS

The combination of

- ◆ a specially designed three-arm hybrid head (combination of basket and grabber),
- ◆ a new Guide wire technology with high shaft stiffness and hyperflexible distal end,
- ◆ including a Tuohy-Borst adapter and an insertion aid for easy insertion into the endoscope's working channel and as splash guard,
- ◆ a wire design designed for durability and functionality,
- ◆ an ergonomic push handle for fatigue-free working and
- ◆ excellent irrigation for best visibility enables the reliable and fast removal of even the smallest concretions from any position

Size (FR/CH)	Length (cm)	Opening Width Ø (mm)	Order Code
1.5	120	8	48012553
1.5	120	11	48012554
2.0	120	8	48012555
2.0	120	11	48012556

Distributed by UROTECH.