



UROTECH®
A MEDI-GLOBE COMPANY




**SUPRAPUBIC BLADDER
DRAINAGE**

SUPRAPUBIC BLADDER DRAINAGE

SUPRA-FIX Puncture- and Exchange Set With pigtail, with funnel

SET COMPONENTS



Pigtail Catheter

- ◆ Polyurethane (PUR)
- ◆ Flexible pigtail to avoid irritation
- ◆ Funnel for attachment of the urinary drainage bag
- ◆ Markings for exact positioning



Split Cannula

- ◆ Medical stainless steel
- ◆ Easy and safe splitting

Guide wire (only in Exchange Set)

- ◆ PTFE-coated

SUPRAPUBIC BLADDER DRAINAGE


Puncture Set / Catheter					
Size (FR/CH)	Length (cm)	Order Code only Pigtail Catheter	Order Code Puncture Set	Split Cannula Size Ø (mm)	Split Cannula Length (cm)
8.0	65	SU-130865	SU-101040	3.6	12
10	65	SU-131065	SU-101046	4.6	12
12	65	SU-131265	SU-101240	5.6	12

Exchange Set						
Size (FR/CH)	Length (cm)	Order Code Exchange Set	Guide wire Length (cm)	Guide wire Diameter	Guide wire Tip	Guide wire Stiffness
8.0	65	SU-131040	100	0.038"	straight	standard
10	65	SU-131046	100	0.038"	straight	standard
12	65	SU-131240	100	0.038"	straight	standard

SUPRAPUBIC BLADDER DRAINAGE

SUPRA-FIX TUR TUR-Puncture Set With pigtail, with funnel

SET COMPONENTS



Pigtail Catheter

- ◆ Polyurethane (PUR)
- ◆ Open or closed tip
- ◆ Very suitable for drainage and irrigation during TUR procedure
- ◆ Flexible pigtail to avoid irritation
- ◆ Optimal urine drainage due to drainage holes in pigtail
- ◆ Funnel for attachment of the urinary drainage bag
- ◆ Markings for exact positioning



Split Cannula

- ◆ Medical stainless steel
- ◆ Easy and safe splitting

SUPRAPUBIC BLADDER DRAINAGE

Open tip

Size (FR/CH)	Length (cm)	Order Code	Split Cannula Size Ø (mm)	Split Cannula Length (cm)
14	65	SU-121640	5.6	12


Closed tip

Size (FR/CH)	Length (cm)	Order Code	Split Cannula Size Ø (mm)	Split Cannula Length (cm)
14	65	SU-121540	5.6	12

SUPRAPUBIC BLADDER DRAINAGE

SUPRA-FIX BK Puncture-, Exchange- and Dilation Exchange Set 100% silicone, with pigtail, with integrated balloon

SET COMPONENTS



Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Flexible pigtail to avoid irritation
- ◆ Markings for exact positioning



Split Cannula

- ◆ Medical stainless steel
- ◆ Easy and safe splitting

Guide wire (only in Exchange Set)

- ◆ PTFE-coated

SUPRAPUBIC BLADDER DRAINAGE

Puncture Set / Catheter

Size (FR/CH)	Length (cm)	Balloon (ml)	Order Code only Balloon Catheter	Order Code Puncture Set	Split Cannula Size Ø (mm)	Split Cannula Length (cm)
12	40	05	SU-251205	SU-121240	5.6	12
14	40	05	SU-251405	SU-121440	5.6	12
16	40	10	SU-251605			
18	40	10	SU-251805			

Exchange Set

Size (FR/CH)	Length (cm)	Balloon (ml)	Order Code Exchange Set	Guide wire Length (cm)	Guide wire Diameter	Guide wire Tip	Guide wire Stiffness
12	40	05	SU-124540	100	0.038"	straight	standard
14	40	05	SU-124640	100	0.038"	straight	standard
16	40	10	SU-124740	100	0.038"	straight	standard
18	40	10	SU-124840	100	0.038"	straight	standard


Dilation Exchange Set

Size (FR/CH)	Length (cm)	Balloon (ml)	Order Code Dilation Exchange Set	Dilator (FR/CH)	Guide wire Length (cm)	Guide wire Diameter	Guide wire Tip	Guide wire Stiffness
12	40	05	SU-124541	14	100	0.038"	straight	standard
14	40	05	SU-124641	16	100	0.038"	straight	standard
16	40	10	SU-124741	18	100	0.038"	straight	standard

SUPRAPUBIC BLADDER DRAINAGE


SUPRA-FIX-INTEGRALE Puncture-, Exchange- and Dilation Exchange Set 100% silicone, with integrated balloon

SET COMPONENTS



Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Markings for exact positioning



Split Cannula

- ◆ Medical stainless steel
- ◆ Easy and safe splitting

Guide wire (only in Exchange Set)

- ◆ PTFE-coated

SUPRAPUBIC BLADDER DRAINAGE

Puncture Set / Catheter						
Size (FR/CH)	Length (cm)	Balloon (ml)	Order Code only Balloon Catheter	Order Code Puncture Set	Split Cannula Size Ø (mm)	Split Cannula Length (cm)
08	30	1.5	SU-270825	SU-110826	3.6	08
10	30	3.0	SU-271025	SU-111026	4.6	08
10	40	3.0	SU-271005	SU-111030	4.6	12
12	30	5.0	SU-271225	SU-111226	4.6	08
12	40	5.0	SU-271205	SU-111230	4.6	12
14	30	5.0	SU-271425			
14	40	5.0	SU-271405	SU-111430	5.6	12
16	30	5-10	SU-271625			
16	40	5-10	SU-271605			
18	30	5-10	SU-271825			
18	40	5-10	SU-271805			
20	30	5-10	SU-272025			
20	40	5-10	SU-272005			
22	30	5-10	SU-272225			
22	40	5-10	SU-272205			
24	40	5-10	SU-272405			

Exchange Set							
Size (FR/CH)	Length (cm)	Balloon (ml)	Order Code Exchange Set	Guide wire Length (cm)	Guide wire Diameter	Guide wire Tip	Guide wire Stiffness
08	30	1.5	SU-140826	100	0.025"	straight	standard
10	30	3.0	SU-141026	100	0.032"	straight	standard
10	40	3.0	SU-141030	100	0.032"	straight	standard
12	30	5.0	SU-141226	100	0.032"	straight	standard
12	40	5.0	SU-141230	100	0.032"	straight	standard
14	30	5.0	SU-141426	100	0.038"	straight	standard
14	40	5.0	SU-141430	100	0.038"	straight	standard
16	30	5-10	SU-141626	100	0.038"	straight	standard
16	40	5-10	SU-141630	100	0.038"	straight	standard
18	30	5-10	SU-141826	100	0.038"	straight	standard
18	40	5-10	SU-141830	100	0.038"	straight	standard
20	40	5-10	SU-142030	100	0.038"	straight	standard


Dilation Exchange Set								
Size (FR/CH)	Length (cm)	Balloon (ml)	Order Code Dilation Exchange Set	Dilator (FR/CH)	Guide wire Length (cm)	Guide wire Diameter	Guide wire Tip	Guide wire Stiffness
10	40	3.0	SU-141031	12	100	0.032"	straight	standard
12	40	5.0	SU-141231	14	100	0.032"	straight	standard
14	40	5.0	SU-141431	16	100	0.038"	straight	standard
16	40	5-10	SU-141631	18	100	0.038"	straight	standard

SUPRAPUBIC BLADDER DRAINAGE

SUPRA-FIX Suprapubic Balloon Catheter

100% silicone, with hollow stylet, with integrated balloon

SET COMPONENTS



Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Reinforced tip
- ◆ Markings for exact positioning



Hollow Stylet

- ◆ To stabilise the balloon catheter

Size (FR/CH)	Length (cm)	Balloon (ml)	Order Code	Hollow Stylet Max. Guide wire Diameter
12	40	5	SU-241205	0.035"
14	40	5	SU-241405	0.035"
16	40	10	SU-241605	0.035"
18	40	10	SU-241805	0.035"