



**UROTECH<sup>®</sup>**  
A MEDI-GLOBE COMPANY




**URETERAL CATHETERS**

Ureteral stents  
nephrostomy  
Suprapubic Bladder drainage  
Accessories  
Ureteral Catheters  
Foley catheters  
Stone baskets  
Special Products

## URETERAL CATHETERS

### URETERAL CATHETERS For the retrograde urography



#### Nelaton tip

Size (FR/CH)	Length (cm)	Order Code closed tip	Order Code open tip	Max. Guide wire Diameter
03	70	UK-500370	UK-600370	0.025"
04	70	UK-500470	UK-600470	0.035"
05	70	UK-500570	UK-600570	0.035"
06	70	UK-500670	UK-600670	0.035"
07	70	UK-500770	UK-600770	0.035"

#### Tiemann tip

Size (FR/CH)	Length (cm)	Order Code closed tip	Order Code open tip	Max. Guide wire Diameter
03	70	UK-510370	UK-540370	0.025"
04	70	UK-510470	UK-540470	0.035"
05	70	UK-510570	UK-540570	0.035"
06	70	UK-510670	UK-540670	0.035"
07	70	UK-510770	UK-540770	0.035"

#### Ureteral catheter by Prof. Walz

Size (FR/CH)	Length (cm)	Tip by Prof. Walz (FR/CH)	Order Code open tip both sides	Max. Guide wire Diameter
04	70	07	UK-630470	0.035"