


Diagonal	
Number	Length (L)
30-41006	60 cm
30-41007	70 cm
30-41008	80 cm
30-41009	90 cm
30-41010	100 cm
30-41011	110 cm
30-41012	120 cm

Tegn. sj		15-11-2012
DS/ISO 2768 m		
Vægt: N/A		N/A
Tegning nr.:		
DiagonalElectric 60-120 cm Reservelede		
ROPOX		
Better ways to better days Ringstedslede 221, 4700 Nestved Tlf.: 53753500 Fax.: 53753539		
Side 1 af 2		

ITEM	QTY	TITLE	3	4	5	6
13	1	Sikkerhedsstop	30-41006	30-41007	30-41009	30-41011
14	1	Topplade	30-41606	30-41607	30-41609	30-41611
15	1	Bundplade	30-40999-104	30-40999-135	30-40999-139	30-40999-143
			30-40999-144	30-40999-125	30-40999-127	30-40999-129
				30-40999-126	30-40999-128	30-40999-115




PART NUMBER	QTY	TITLE	L	B	ITEM
96000592	1	Netkabel			17
30-40999-160	2	Rør til topbeslag			16
98002035	1	Modstand 2,2 K ohm			12
96000593	1	Motor kabel			11
98002008	1	Ledning til sikkerhedsstop			10
98002082	1	Adapterstik til sikkerhedsstop			9
96000598	1	Adapterstik til controlboks			8
30-67840-7	1	Bejening			7
30-40999-237	1	Motor med spindel			6
97701-832	1	plade til controlboks			5
96000591	1	Controlboks			4
30-40601-144	1	Plade til ledninger			3
97701-732	1	controlboksplade holder			2
30-40902-000	1	Diagonal			1

Tegn. nr.		15-11-2012
DS/ISO		2768 m
Vægt:		N/A
Tegning nr.:		

DiagonalElectric 60-120 cm
Reserve dele



Better ways to better days
Ringsledgade 221, 4700 Næstved
Tlf. 55790500 Fax. 55790550



Number	Length (L)
30-41013	130 cm
30-41014	140 cm
30-41015	150 cm
30-41016	160 cm
30-41017	170 cm
30-41018	180 cm

Tegn. nr. 19-11-2012	
DS/ISO 2768 m	
Vægt: N/A	N/A
Tegning nr.:	
DiagonalElectric 130-180 cm Reserve dele	
<p>ROPOX</p> <p>Better ways to better days Blåstredgade 221, 4700 Næstved Tlf.: 55750500 Fax.: 55750550</p>	
Side 1 af 2	

